

---

ReadyAgarose™ Gels  
Instruction Manual

Catalog #  
161-3000

**Bio-Rad Laboratories, Inc.**

2000 Alfred Nobel Dr., Hercules, CA 94547 USA  
510-741-1000

4110108 Rev C



# Table of Contents

<b>Section 1 General Information .....</b>	<b>1</b>
1.1 ReadyAgarose Gel Specifications .....	1
1.2 ReadyAgarose Comb Configurations .....	2
<b>Section 2 Basic Operation .....</b>	<b>3</b>
2.1 Setting Up and Running ReadyAgarose Gels .....	3
2.2 Removing the Auto-Lock Tabs .....	6
<b>Section 3 Electrophoresis Information.....</b>	<b>7</b>
3.1 Electrophoresis Buffers .....	7
3.2 Relative Sample Migration Rates.....	8
3.3 Loading Dye Migration Characteristics .....	8
3.4 Submerged Horizontal Electrophoresis Subcell Compatibility .....	9
<b>Section 4 Troubleshooting .....</b>	<b>10</b>
<b>Section 5 Ordering Information.....</b>	<b>11</b>

# Section 1 General Information

## 1.1 ReadyAgarose Gel Specifications

Gel material	1% gels: Certified™ molecular biology agarose 3% gels: Certified low range ultra agarose	
Gel dimensions (W x L)	Mini gel: 6.6 cm x 10 cm; wide mini gel: 15.1 cm x 10 cm; 96 Plus gel: 15.1 cm x 10 cm	
Gel thickness	5.5 mm	
Optimal resolution range	<b>Mini Gels and Wide</b>	
<b>Gel percentage</b>	<b>Mini Gels</b>	<b>96 Plus Gels</b>
1%	200–10,000 bp	500–10,000 bp
3%	20–2,000 bp	20–2,000 bp
Tray dimensions	Mini tray: 7.2 x 10 cm; wide mini tray: 15.6 x 10 cm; 96 Plus tray: 15.6 x 10 cm	
Tray material	UV-transparent acrylic	
Packaging tray lid material	PETG	
Storage conditions	Store label side up at 4–22°C; DO NOT FREEZE.	

## 1.2 ReadyAgarose Comb Configurations


Gel Type	Comb	Load Volume	Comments
Mini	8-well	20 $\mu$ l	
Mini	12-well	10 $\mu$ l	
Mini	2 x 8-well	20 $\mu$ l	
Wide mini	20-well	20 $\mu$ l	
Wide mini	32-well	15 $\mu$ l	Multichannel pipet compatible, 12 $\mu$ l
Wide mini	2 x 32-well	15 $\mu$ l	Multichannel pipet compatible, 12 $\mu$ l
96 Plus	4 x 26-well	15 $\mu$ l	Multichannel pipet compatible, 12 $\mu$ l

2

## Section 2 Basic Operation

### 2.1 Setting Up and Running ReadyAgarose Gels

1. Each ReadyAgarose gel is packaged individually in a foil bag. Open the bag by pulling at the tear notch. Carefully remove the gel from the foil bag. Remove the protective gel tray cover.



3


- 2.** Place the gel in the recommended Bio-Rad Sub-Cell® GT cell. The two auto-lock tabs on the tray fit into the grooves on the sides of the Bio-Rad Sub-Cell GT, which correctly positions the tray. Make sure the wells are oriented toward the cathode (black negative electrode). ReadyAgarose gels are designed to work with most subcells. See section 2.2 if you need to remove the auto-lock tabs and section 3.4 for a list of compatible subcells.
- 3.** Fill the Sub-Cell GT with electrophoresis buffer so that 2–3 mm of buffer is above the gel. For best results, do not fill buffer above the sides of the tray. For ReadyAgarose gels containing ethidium bromide, additional ethidium bromide can be added to the running buffer at a concentration of 0.5 µg/ml to avoid faint bands near the bottom of the gel. To ensure 2–3 mm of buffer over the gel when using a Bio-Rad Sub-Cell GT, use 210 ml buffer with the Mini-Sub® cell GT and 550 ml with the wide Mini-Sub cell GT.
- 4.** Load samples and standards. For the ReadyAgarose 96 Plus gels, load the DNA standard or marker in the outer two wells (labeled “M” in the top row).
- 5.** Place the lid on the Sub-Cell GT, aligning the color-coded banana plugs and jacks. Set voltage conditions and start the run. To obtain the best results with

ReadyAgarose gels, a voltage gradient of 3–6.5 V/cm (measured as the distance between the electrodes) should be used (A setting of 75 V for the Mini-Sub cell GT and wide Mini-Sub cell GT creates a voltage gradient of 5 V/cm).

- 6.** After electrophoresis is complete, turn off the power supply and disconnect the electrical leads. Remove the gel from the Sub-Cell.
- 7.** For ReadyAgarose gels without ethidium bromide, a gel stain is required to visualize the nucleic acids in the gel. We suggest using the following method for staining with ethidium bromide (for staining with other fluorescent stains or silver, follow the manufacturer’s instructions). Transfer the agarose gel to a staining dish with 0.5 µg/ml ethidium bromide. Stain the gel for 15–30 min. Remove excess stain from the gel by placing it in water for 10–30 min.
- 8.** ReadyAgarose gels can be viewed by placing the UV-transparent tray containing the agarose gel directly on a UV transilluminator. Alternatively, the agarose gel can be removed from the tray before UV visualization for maximum sensitivity.

## 2.2 Removing the Auto-Lock Tabs

1. To use ReadyAgarose gels with non-Bio-Rad Subcell GT cells, one or both of the auto-lock tabs may need to be removed for the tray to fit into these cells.
2. To remove the auto-lock tabs, place the “gel” (tray with gel) at approximately a 45° angle, with one length of the tray touching the bench. Place a finger over the auto-lock tab and press down to break the tab off of the tray as shown in the figure above.



6

## Section 3 Electrophoresis Information

### 3.1 Electrophoresis Buffers

DNA agarose electrophoresis is usually conducted with either Tris-acetate-EDTA (TAE) or Tris-borate-EDTA (TBE) buffer. While TAE provides faster electrophoretic migration of linear DNA, TBE buffers have a stronger buffering capacity for longer or higher voltage electrophoresis runs.

#### Buffer

1x TAE buffer  
1x TBE buffer  
5x Nucleic acid sample buffer

#### Formulation

40 mM Tris, pH 7.6, 20 mM acetic acid, 1 mM EDTA  
89 mM Tris, pH 7.6, 89 mM boric acid, 2 mM EDTA  
50 mM Tris, 25% glycerol, 5 mM EDTA,  
0.2% Bromophenol Blue, 0.2% Xylene Cyanole FF

7

### 3.2 Relative Sample Migration Rates

<b>Voltage</b>	<b>Voltage Gradient</b>	<b>Bromophenol Blue Migration Rate</b>
30 V	2 V/cm	1.2 cm/hr
45 V	3 V/cm	1.9 cm/hr
60 V	4 V/cm	2.8 cm/hr
75 V	5 V/cm	3.8 cm/hr
90 V	6 V/cm	4.9 cm/hr
105 V	7 V/cm	5.5 cm/hr

### 3.3 Loading Dye Migration Characteristics\*

<b>Agarose Concentration</b>	<b>Xylene Cyanole</b>	<b>Bromophenol Blue</b>
1%	6,000 bp	700 bp
3%	300 bp	50 bp

\*Approximate

### 3.4 Submerged Horizontal Electrophoresis Subcell Compatibility

ReadyAgarose gels are designed to fit most subcells. ReadyAgarose mini gels will fit subcells that are 7.1 cm wide (inside dimension). ReadyAgarose wide mini gels will fit subcells that are 15.6 cm wide (inside dimension). ReadyAgarose gels are compatible with the subcells listed below.

<b>Company</b>	<b>ReadyAgarose Mini Gel</b>	<b>ReadyAgarose Wide Mini Gel and ReadyAgarose 96 Plus Gel</b>
Bio-Rad	Mini Sub-Cell GT	Wide Mini-Sub cell GT
Owl Separations	EasyCast Minigel B1A, B1, B2, Buffer Puffer	Centipede wide format, Gator large format A2, Gator large format A3
C.B.S. Scientific	MGU-202T, MGU-252T, MGU-402T	Horizontal gel system
Amersham	Hoefer HE33 mini, Hoefer HE99X	GNA 200
Stratagene	Joule Box Horizontal elec. app.	Horizontal electrophoresis apparatus

## Section 4 Troubleshooting

<b>Problem</b>	<b>Possible Cause</b>	<b>Solution</b>
Curved bands, smiles	Sample overload Glycerol concentration is too high	Reduce total sample volume Reduce final glycerol concentration to ~ 5%
Band smearing	Excessive voltage and heating	1. Reduce voltage 2. Check buffer composition
	Sample spilled out of well	1. Apply sample carefully 2. Ensure correct loading dye concentration is used
Slanted lanes	Salt concentration in sample too high	Reduce salt concentration
	Subcell not level	Level subcell or place on leveling table
Broad distorted bands	Sample concentration too high	Reduce total concentration of nucleic acid in the sample

**10**

## Section 5 Ordering Information

### ReadyAgarose Mini Gels

<b>TBE Mini Gels</b>	<b>8-Well</b>	<b>12-Well</b>
1% plus EtBr*	161-3004	161-3010
3% plus EtBr*	161-3006	161-3012

<b>TAE Mini Gels</b>	<b>8-Well</b>	<b>12-Well</b>
1%	161-3015	
1% plus EtBr*	161-3016	161-3022
3%	161-3017	
3% plus EtBr*	161-3018	161-3024
1%	161-3057 (2 X 8-well)	

### ReadyAgarose 96 Plus Gels

<b>96 Plus Gels</b>	<b>TBE</b>
1% plus EtBr*	161-3060
3% plus EtBr*	161-3062

### ReadyAgarose Wide Mini Gels

<b>TBE Wide Mini Gels</b>	<b>20-Well</b>	<b>32-Well</b>	<b>2 x 32-Well</b>
1% plus EtBr*	161-3028	161-3034	161-3038
3% plus EtBr*	161-3030	161-3036	161-3040

<b>TAE Wide Mini Gels</b>	<b>20-Well</b>	<b>32-Well</b>	<b>2 x 32-Well</b>
1%			
1% plus EtBr*	161-3044	161-3050	161-3054
3%			
3% plus EtBr*	161-3046	161-3052	161-3056

### TAE

161-3063
161-3065

\*Ethidium bromide

**11**

**ReadySub-Cell™ Systems**

- 170-4487 Mini ReadySub-Cell GT cell, for ReadyAgarose gels, includes subcell unit only  
170-4489 Wide Mini ReadySub-Cell GT cell, for ReadyAgarose gels, includes subcell unit only

**Electrophoresis Buffers and Stains**

- 161-0743 50x Tris/Acetic Acid/EDTA, 1 L  
161-0773 50x Tris/Acetic Acid/EDTA, 5 L cube  
161-0733 10x Tris/Boric Acid/EDTA, 1 L  
161-0770 10x Tris/Boric Acid/EDTA, 5 L cube  
161-0774 20x SSC, 1 L  
161-0775 20x SSC, 5 L cube  
161-0433 Ethidium Bromide Solution, 10 ml, 10 mg/ml  
161-0430 Ethidium Bromide Tablets, 10 x 10 mg

**Nucleic Acid Standards**

- 170-8351 EZ Load 20 bp Molecular Ruler  
170-8352 EZ Load 100 bp Molecular Ruler  
170-8353 EZ Load 100 bp PCR Molecular Ruler  
170-8200 AmpliSize® Molecular Ruler  
170-8354 EZ Load 500 bp Molecular Ruler  
170-8355 EZ Load 1.0 kb Molecular Ruler  
170-8205 EZ Load 2.5 bp Molecular Ruler

Buffer Puffer, Centipede, Easy Cast and Gator are trademarks of Owl Separation Systems.

Hoefer is a trademark of GE Healthcare.

Joule Box is a trademark of Stratagene.