

# Bio-Plex™ Pro Human Cytokine, Chemokine, and Growth Factor Assays

## Quick Guide

For Use with	Instruction Manual #
Bio-Plex Pro Human Cytokine, Chemokine, and Growth Factor Assays	10000111560

This guide can be used to prepare and run a full 1 x 96-well assay plate. New users can go to [bio-rad.com/HCS](http://bio-rad.com/HCS) and download the manual, which includes detailed instructions and a list of kit components.

**IMPORTANT!** Pay close attention to **vortexing**, **shaking**, and **incubation** instructions. Deviation from the protocol may result in low assay signal and assay variability.

## Initial Preparation

1. Plan the plate layout.
2. Start up/warm up the Bio-Plex Multiplex Immunoassay System (**30 min**).
  - Bring diluents, including wash buffer, assay buffer, standard diluent HB, detection antibody diluent HB, and sample diluent HB, to room temperature (RT). Keep the other items on ice until needed
    - Mix by inversion to ensure all salts are in solution
    - Prepare 1x wash buffer: dilute **1 part** 10x wash buffer (60 ml) with **9 parts** distilled water (540 ml)
  - Begin to thaw frozen samples
3. Prepare the sample dilution according to the guidelines provided in the following table. It is important to centrifuge serum or plasma samples at **1,000 x g** for **15 min** at **4°C** to remove particulates from all samples prior to use.


Sample Type	Recommended Sample Dilution	Diluent
Serum and plasma	1:4	Sample diluent HB
Culture media and fluids	User defined	Diluent + 0.5% bovine serum albumin (w/v)

**Note:** ICAM-1 and VCAM-1 require higher dilution for serum and plasma (100-fold recommended). Refer to the Bio-Plex Pro Human Cytokine, Chemokine, and Growth Factor Assays Instruction Manual (10000111560) for detailed sample preparation recommendations.

## Bio-Plex Pro Human Cytokine, Chemokine, and Growth Factor Assays Quick Guide

- Prime the wash station for a flat bottom plate or set a vacuum manifold to -1 to -3" Hg for filter plate.
- Calibrate the Bio-Plex System by following the prompts in Bio-Plex Manager Software.
- Reconstitute the standards and control by adding 250  $\mu$ l of standard diluent HB to each. **Vortex** at medium speed for **5 sec** and incubate all vials on ice for precisely **30 min**.
- Prepare a fourfold standard dilution series and blank as shown. **Vortex** at medium speed for **5 sec** between liquid transfers.

**Note:** Standards are at S1 concentration after reconstitution. Controls are ready to use after reconstitution and no dilution is needed. Controls are included with the fixed panel only.



- Vortex** the coupled beads at medium speed for **30 sec** and **dilute to 1x** in Bio-Plex Assay Buffer as shown. Protect from light.

### Premixed Panels

Number of Wells	10x Beads, $\mu$ l	Assay Buffer, $\mu$ l	Total Volume, $\mu$ l
96	575	5,175	5,750

### Singleplex Assays

Number of Wells	Singleplex #1	Singleplex #2	Assay Buffer, $\mu$ l	Total Volume, $\mu$ l
	20x Beads, $\mu$ l	20x Beads, $\mu$ l		
96	288	288	5,184	5,760

**Note:** 20x singleplex beads allow multiplexing up to 20 analytes.

## Running the Assay

- Vortex** the diluted (1x) beads. **Add 50 µl** to each well of the assay plate.
- Wash the plate two times** with **100 µl** Bio-Plex Wash Buffer.
- Vortex** the samples, standards, blank, and controls. **Add 50 µl** to each well.
- Cover the plate with sealing tape. Incubate on shaker at **850 ± 50 rpm** at RT for **30 min**.
- With 10 min left in the incubation, **vortex** the detection antibodies for **5 sec** and quick-spin to collect liquid. **Dilute to 1x** as shown.
- Wash the plate three times** with **100 µl** wash buffer.

### Premixed Panels

Number of Wells	10x Detection Antibodies, µl	10x Detection Antibody Diluent HB, µl	Total Volume, µl
96	300	2,700	3,000

### Singleplex Assays

Number of Wells	Singleplex #1 20x Detection Antibodies, µl	Singleplex #2 20x Detection Antibodies, µl	Detection Antibody Diluent HB, µl	Total Volume, µl
96	150	150	2,700	3,000

**Note:** 20x singleplex beads allow multiplexing up to 20 analytes.

- Vortex** the diluted (1x) detection antibodies. **Add 25 µl** to each well.
- Cover the plate with sealing tape and incubate at **850 ± 50 rpm** for **30 min** at RT. Meanwhile, prepare the Bio-Plex Manager Software protocol; enter standard S1 values and units provided in the assay kit.
- With 10 min left in the incubation, **vortex** 100x streptavidin-phycoerythrin (SA-PE) for **5 sec** and quick-spin to collect liquid. **Dilute to 1x** as shown and protect from light.

Number of Wells	100x SA-PE, µl	Assay Buffer, µl	Total Volume, µl
96	60	5,940	6,000

- Wash the plate three times** with **100 µl** wash buffer.
- Vortex** the diluted (1x) SA-PE. **Add 50 µl** to each well.

12. Cover the plate with sealing tape and incubate at **850 ± 50 rpm** for **10 min** at RT.
13. Wash the plate three times with **100 µl** wash buffer.
14. Resuspend the beads in **125 µl** assay buffer. Cover and shake at **850 ± 50 rpm** for **30 sec**.
15. Remove the sealing tape and **read plate** using the following settings:

Instrument	RP1 (PMT)	DD Gates	Bead Events
Bio-Plex 3D*	Standard	Select MagPlex Beads	50
Bio-Plex 100, 200*	Low	5,000–25,000	50
Luminex MAGPIX	N/A, use default instrument settings		

\* Or similar Luminex System.

16. Controls are included with the fixed panel only. If the control was run, compare the observed concentration against the ranges provided in the assay kit. Ranges apply only when standards and controls are prepared in Bio-Plex Standard Diluent HB.

BIO-RAD and BIO-PLEX are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. Luminex is a trademark of Luminex Corporation. All trademarks used herein are the property of their respective owner.  
© 2023 Bio-Rad Laboratories, Inc.

The Bio-Plex Suspension Array System includes fluorescently labeled microspheres and instrumentation licensed to Bio-Rad Laboratories, Inc. by the Luminex Corporation.

## BIO-RAD

**Bio-Rad  
Laboratories, Inc.**

Life Science  
Group

**Website** bio-rad.com **USA** 1 800 424 6723 **Australia** 61 2 9914 2800  
**Austria** 00 800 00 24 67 23 **Belgium** 00 800 00 24 67 23 **Brazil** 4003 0399  
**Canada** 1 905 364 3435 **China** 86 21 6169 8500 **Czech Republic** 00 800 00 24 67 23  
**Denmark** 00 800 00 24 67 23 **Finland** 00 800 00 24 67 23 **France** 00 800 00 24 67 23  
**Germany** 00 800 00 24 67 23 **Hong Kong** 852 2789 3300 **Hungary** 00 800 00 24 67 23  
**India** 91 124 4029300 **Israel** 0 3 9636050 **Italy** 00 800 00 24 67 23  
**Japan** 81 3 6361 7000 **Korea** 82 080 007 7373 **Luxembourg** 00 800 00 24 67 23  
**Mexico** 52 555 488 7670 **The Netherlands** 00 800 00 24 67 23 **New Zealand** 64 9 415 2280  
**Norway** 00 800 00 24 67 23 **Poland** 00 800 00 24 67 23 **Portugal** 00 800 00 24 67 23  
**Russian Federation** 00 800 00 24 67 23 **Singapore** 65 6415 3188  
**South Africa** 00 800 00 24 67 23 **Spain** 00 800 00 24 67 23 **Sweden** 00 800 00 24 67 23  
**Switzerland** 00 800 00 24 67 23 **Taiwan** 886 2 2578 7189 **Thailand** 66 2 651 8311  
**United Arab Emirates** 36 1 459 6150 **United Kingdom** 00 800 00 24 67 23

