

Dark Probe Mutation Detection Assays

More Targets. Less Sample. Higher Sensitivity.

Dark Probe assays decouple wild-type suppression from fluorescence, freeing detection channels on Bio-Rad Droplet Digital™ PCR (ddPCR™) Systems.* Detect more mutations per reaction, and include reference genes for accurate VAF calculation, without increasing sample input or sacrificing sensitivity.

How It Works

Replace the wild-type fluorescent probe with a quencher-only dark probe to maintain specificity while freeing each detection channel for an additional mutation target or reference gene. The freed channels can be used to increase assay multiplexing, incorporate reference genes for VAF normalization, or both, all from the same sample input.

Standard vs. Dark Probe Approach

Standard mutation assay	Dark probe approach
• 2 fluorescent probes (wild-type and mutant) • 2 channels per mutation • Maximum 3 mutations per well	• Non-fluorescent WT probes • 1 channel per mutation target • Up to 4 mutations and 2–3 reference genes per well (system-dependent) • Include VAF reference genes (e.g., <i>RPP30</i>, <i>AGO1</i>) without a separate reaction

Get More Targets From Every Well

Instrument	Standard Mutation Detection Assays	Dark Probe Mutation Detection Assays
QX200™ ddPCR System†	1 mutation and 1 WT per well	2 mutations or 1 mutation and 1 reference gene per well
QX600 ddPCR System	3 mutations per well	4 mutations and 2 reference genes per well
QX700™ ddPCR System*	3 mutations per well	4 mutations and 3 reference genes per well

Reference genes (e.g., *RPP30*, *AGO1*) can be included in freed channels for VAF/FA% normalization. No additional reaction required.

Key Benefits



Higher multiplex density:

detect up to four mutations plus up to three reference genes in a single well (system-dependent), without increasing sample input.



Lower cost per target:

consolidate panels into one well to reduce wells, consumables, reagent volume, and preserve precious sample.



Flexible panel design:

no required fluorophore pairings (FAM/HEX, Cy5/Cy5.5, ROX/ATTO). Freely assign channels to mutations or reference genes based on your panel design. Order through the Bio-Rad Mutation Detection Design Engine.

Available Formats

Channels: FAM, HEX, Cy5, Cy5.5 Reaction sizes: 200 x 20 µL, 1,000 x 20 µL, 2,500 x 20 µL

* For Research Use Only. Not for use in diagnostic procedures.

† For Research Use Only. In vitro diagnostic (IVD) option is available.



Visit [bio-rad.com/MutationDetection](https://www.bio-rad.com/MutationDetection) for more information.

Common Applications

- Oncology research:** mutation profiling panels for *EGFR*, *KRAS*, *BRAF*, and other variants across tissue, FFPE, and ctDNA samples.
- MRD and liquid biopsy:** track rare variants and monitor therapy responses at low VAF in ctDNA.
- Translational research:** confirm NGS calls and quantify SNVs in a single reaction.

Performance At a Glance

- Sensitivity:** 0.05–0.1%+ VAF at 100 ng input per well; lower VAF achievable with increased input or multiple wells.
- Multiplexing:** six targets: four mutations and two reference genes can be resolved in a single well (Figure 1).

+ Achievable sensitivity is driven by sample input, target abundance, and assay design.




Fig. 1. Multiplex detection of four mutations and two reference genes in a single QX600 ddPCR System well. Six targets resolved across three channel pairs in a single QX600 ddPCR System well: four mutations in panels **A** and **B** (with dark probes suppressing wild-type signal), and two reference genes co-amplifying in panel **C** for normalization.

Ordering Information

Channel	Description	200 x 20 µL Reactions	1,000 x 20 µL Reactions	2,500 x 20 µL Reactions
FAM	ddPCR MUT FAM+Dark Assay	12016634	12016618	12016619
HEX	ddPCR MUT HEX+Dark Assay	12016607	12016586	12016617
Cy5	ddPCR MUT Cy5+Dark Assay	12016608	12016638	12016624
Cy5.5	ddPCR MUT Cy5.5+Dark Assay	12016635	12016636	12016637

ROX and ATTO 590 assays: coming soon. All Dark Probe assay sequences are designable via the Bio-Rad Mutation Detection Design Engine.

Visit [bio-rad.com/MutationDetection](https://www.bio-rad.com/MutationDetection) for more information.

Bio-Rad, QX200, QX600, QX700, Droplet Digital, Droplet Digital PCR, and ddPCR are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner. © 2026 Bio-Rad Laboratories, Inc.

Purchase of Digital PCR and/or Single-Cell NGS Sample Preparation products (the "Products") from Bio-Rad Laboratories is subject to Bio-Rad Laboratories, Inc. Standard Terms and Conditions of Sale, which can be accessed at <https://www.bio-rad.com/en-us/terms-conditions>. Unless we expressly state otherwise in additional Terms and Conditions, no rights are granted for you to distribute or resell the Products. Unless we expressly state otherwise in additional Terms and Conditions, no rights are granted for the development or commercialization of diagnostic assays for use with the Products without a license from Bio-Rad. It is the user's obligation to obtain a commercial license from Bio-Rad for (i) all commercial uses (not just diagnostic uses) and (ii) sale of assays for use on Bio-Rad's dPCR and ddSEQ instruments. The Products and/or their use are covered by U.S. and foreign patents and/or pending patent applications owned by or under license to Bio-Rad Laboratories, Inc. See <https://www.bio-rad.com/en-us/trademarks>.



Bio-Rad
Laboratories, Inc.

Life Science Group	Website	bio-rad.com	USA	1 800 424 6723	Australia	61 2 9914 2800	Austria	00 800 00 24 67 23	Belgium	00 800 00 24 67 23	Brazil	55 11 3065 7550
	Canada	1 800 361 1808	China	86 21 6169 8500	Czech Republic	00 800 00 24 67 23	Denmark	00 800 00 24 67 23	Finland	00 800 00 24 67 23	France	00 800 00 24 67 23
	Germany	00 800 00 24 67 23	Greece	30 210 7774396	Hong Kong	852 2789 3300	Hungary	00 800 00 24 67 23	India	91 124 4029300	Israel	000 800 00 24 67 23
	Italy	00 800 00 24 67 23	Japan	81 3 6361 7000	Korea	82 080 007 7373	Luxembourg	00 800 00 24 67 23	Mexico	52 55 5488 7670	The Netherlands	00 800 00 24 67 23
	New Zealand	64 9 415 2280	Norway	00 800 00 24 67 23	Poland	00 800 00 24 67 23	Portugal	00 800 00 24 67 23	Russian Federation	7 495 721 14 04	Singapore	65 6415-3170
	South Africa	27 21 531 7504	Spain	00 800 00 24 67 23	Sweden	00 800 00 24 67 23	Switzerland	00 800 00 24 67 23	Taiwan	886 2 2578 7189	Thailand	662 651 8311
	United Arab Emirates	971 4 818 7300	United Kingdom	00 800 00 24 67 23								