PIF Supplement

Catalog #	Description
12013322	PIF Supplement, 2 g

For laboratory use only.

Intended Use

The PIF Supplement is used with buffered peptone water (BPW) in food microbiology for the enrichment of *Enterobacteriaceae*, particularly *Salmonella* spp. and *Cronobacter* spp. It is adapted for matrices with a high level of background flora, such as lactic acid bacteria, bifidobacteria, and other gram-positive bacteria. It is especially suitable for the detection of *Salmonella* spp. and *Cronobacter* spp. in powdered infant formula with and without probiotics, ingredients for infant formula, and environmental samples of production areas.

Principle

The PIF Supplement is a special mix of components, which allows the inhibition of background flora and, when associated with Bio-Rad's BPW, optimized recovery of *Enterobacteriaceae*, in particular *Salmonella* and *Cronobacter*. The PIF Supplement may be used for the enrichment step in the following kit protocols:

- iQ-Check *Cronobacter* PCR Detection Kit
- iQ-Check *Enterobacteriaceae* PCR Detection Kit
- iQ-Check *Salmonella* II PCR Detection Kit
- RAPID'*Salmonella* Agar
- RAPID'*Sakazakii* Agar

Theoretical Composition

The PIF Supplement is supplied as a specific colored mix in 2 g amounts.

Shelf Life and Storage

Store dehydrated at 2–8°C in a carefully sealed package in a dry and dark place. The PIF Supplement is stable for 6 weeks at 2–8°C in a dark place once reconstituted.

Required Materials Not Supplied

This is a non-exhaustive list.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer

Supplies

- Buffered Peptone Water (catalog #3554179, 225 ml x 6 bottles; 3564684, dehydrated 500 g; 3564686, dehydrated, 5 kg; 3555795, 3 L x 4 bags; 3555790, 5 L x 2 bags)
- iQ-Check *Cronobacter* PCR Detection Kit (catalog #3578137)
- iQ-Check *Enterobacteriaceae* PCR Detection Kit (catalog #12003068)
- iQ-Check *Salmonella* II PCR Detection Kit (catalog #3578123)
- RAPID'*Sakazakii* Agar (catalog #3563971, 90 mm x 20 plates; 3564976, dehydrated, 500 g)

- RAPID' *Salmonella* Agar (catalog #3563961, 90 mm x 20 plates; 3563963, 90 mm x 120 plates 3564705, dehydrated, 500 g)
- Sterile distilled water (catalog #3554154, 9 ml x 25 tubes; 3554283, 18 ml x 60 tubes)
- Sterile Pasteur pipets
- Sterile pipets
- Sterile weigh bags

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria
- Product development and validation of the PIF Supplement have been performed with Bio-Rad's BPW for optimal performance of the method
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- For SDS product safety information and certificate of analysis, visit bio-rad.com

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

PIF Supplement Preparation

- Dissolve 2 g of PIF Supplement in 40 ml of 1:1 ethanol/ sterile distilled water solution directly in the plastic vial Mix, heating if necessary, until a homogeneous suspension is obtained
- Homogenize by vortexing
- Leave on the bench for 5 min at room temperature
- Always mix the reconstituted solution before use
- After reconstitution, the solution can be stored at 2–8°C in a dark place for 6 weeks

Enrichment Protocol

Sample Size, g	Volume of BPW at 1:4 Dilution, ml	Volume of PIF Supplement, μ l
10	30	4
25	75	10
100	300	40
300	900	120
375	1,125	150

- Add n g of samples in $3 \times n$ ml of prewarmed BPW (see table)
- Add $0.1 \times v$ ml of PIF solution to v L of matrix/enrichment broth suspension
- Homogenize
- Incubate without shaking at $37 \pm 1^\circ\text{C}$ for 18–24 hr for RAPID methods or for 18–26 hr for iQ-Check methods
- After the enrichment step, the broth can be stored at 2–8°C for 72 hr before performing the screening step

Screening Protocol

Refer to the appropriate user guide for the protocol of each screening method.

Organism	Screening Method	Scope	Validation
<i>Salmonella</i> spp.	iQ-Check <i>Salmonella</i> II RAPID' <i>Salmonella</i>	Infant formula and cereals with and without probiotics and ingredients for infant formula (50–375 g samples) with 1:4 dilution	NF Validation BRD 07/06-07/04 NF Validation BRD 07/11-12/05
<i>Cronobacter</i> spp.	iQ-Check <i>Cronobacter</i> RAPID' <i>Sakazakii</i>		NF Validation BRD 07/23-01/13 NF Validation BRD 07/22-05/12
<i>Enterobacteriaceae</i>	iQ-Check <i>Enterobacteriaceae</i>	Milk powder, infant formula with and without probiotics (10 g and 375 g), and stainless steel (4 x 4 in.)	AOAC 082003

Revision History

Release date	Document number	Change
April 2020	10000126749 Ver A	New document
June 2020	10000126749 Ver B	- Revision history table added - Translations (French, German, Spanish, Italian, and Portuguese)
June 2021	10000126749 Ver C	Update of storage time in Protocol section
January 2022	10000126749 Ver D	Translation (Simplified Chinese)

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

iQ-CHECK is a trademark of Bio-Rad Europe GmbH in certain jurisdictions. All trademarks used herein are the property of their respective owner.

All trademarks used herein are the property of their respective owner.

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