Half Fraser Standard Broth

Catalog #	Description
12018707	Half Fraser Standard , ready-to-use, 5 L x 2 bags
12018706	Half Fraser Standard , dehydrated, 500 g
12018705	Half Fraser Standard , dehydrated, 5 kg

For laboratory use only.

Intended Use

Half Fraser Standard broth is used for the enrichment of *Listeria monocytogenes* and other *Listeria* spp. in food products and environmental samples. This broth can be used as primary enrichment broth following the ISO 11290-1 or as the single enrichment broth of many alternative methods for the detection of *Listeria monocytogenes* and other *Listeria* spp.

Principle

This broth contains a particular mix of nutrient to optimize the growth of *Listeria* spp. strains including *Listeria monocytogenes*. Due to the combined action of lithium chloride, acriflavine, and nalidixic acid, this broth inhibits non-*Listeria* flora. The balance between nutritivity and selectivity ensures the good growth of *Listeria monocytogenes* and other *Listeria* spp. in all types of matrices.

Theoretical Composition

Peptones	5 g
Tryptone	5 g
Yeast extract	5 g
Meat extract	5 g
Sodium chloride	20 g
Disodium hydrogen phosphate dihydrate	12 g
Potassium dihydrogen phosphate	1.35 g
Esculin	1 g
Lithium chloride	3 g
Acriflavine hydrochloride	12.5 mg
Nalidixic acid	10 mg
Ferric ammonium citrate (III)	500 mg
Distilled water	1,000 ml
Final pH at 25°C = 7.2 ± 0.2	

Shelf Life and Storage

Store ready-to-use and dehydrated media at 15–25°C in a carefully sealed package in a dry and dark place until the expiration date. Store media prepared from dehydrated at 2–8°C for 2 weeks in a dark place.

Required Materials Not Supplied

This list is not exhaustive.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria.
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Bags of broth must be manipulated under aseptic conditions. After opening, the bags can be stored carefully sealed with a sterile clamp at ambient temperature for 3 days. Bags can be re-opened once. Do not use the volume remaining in the spout, remove it before dispensing the broth
- For SDS product safety information and certificate of analysis, visit bio-rad.com

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

Dehydrated Half Fraser Standard Broth Preparation

- Shake the bottle before use
- Dissolve 57.9 g of Half Fraser Standard Broth in 1 L of sterile distilled water
- Heat gently, agitating frequently, then bring to a boil until a homogeneous suspension is obtained
- Dispense 225 ml per bottle
- Sterilize in autoclave at 121 ± 3 °C for 15 min

ISO Sample Preparation and Enrichment Protocol

- Prepare sample according to the standard method applicable to the product concerned
- Homogenize 25 g of test sample in 225 ml of Half Fraser Standard Broth
- Incubate primary enrichment for 24–26 hr at 30 ± 1 °C
- Inoculate the surface of one AL and one Palcam, Oxford, or RAPID'L *mono* plate
- Refer to ISO 11290-1 for agar media incubation conditions
- In addition, transfer 0.1 ml of primary enrichment to a 10 ml tube of Fraser Broth
- Incubate for 22-26 hr at 37 ± 1 °C
- Inoculate the surface of one AL and one Palcam, Oxford, or RAPID'L *mono* plate
- Refer to ISO 11290-1 for agar media incubation conditions

Chromogenic Media Sample Preparation and Enrichment Protocol

Refer to the user guides for AL (document #10000135677) and RAPID'L *mono* Agar (document #10000127436) for enrichment and detection protocols.

References

Curtis GDW et al. (1989). A selective differential medium for the isolation of *Listeria monocytogenes*. Lett Appl Microbiol 8, 95–98.

Curtis GDW et al. (1989). Selective agents for *Listeria* can inhibit their growth. Lett Appl Microbiol 8, 169–172.

ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 1: Detection method.

ISO 11290-2:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 2: Enumeration method.

Lovett J et al. (1987). *Listeria monocytogenes* in raw milk: Detection, incidence, and pathogenicity. J Food Prot 50, 188–192.

Prentice GA and Neaves P (1988). *Listeria monocytogenes* in Food: Its Significance and Methods for Its Detection. (Brussels: International Dairy Federation), p. 223.

Van Netten P et al. (1988). A selective and diagnostic medium for use in the enumeration of *Listeria* spp. In foods. Int J Food Microbiol 6, 187–198.

Revision History

Release date	Document number	Change
October 2022	10000159301 Ver A	New document

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

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