iQ-Check *Cronobacter* spp. Kit

User Guide

Test for the real-time PCR detection of *Cronobacter* spp. in infant formula, infant cereal, and environmental samples

Catalog #3578137

BIO-RAD

Table of Contents

Section 1. Introduction	1
Section 2. iQ-Check <i>Cronobacter</i> spp. Technology	1
Section 3. Kit Components.....	2
Section 4. Shelf Life and Storage	3
Section 5. Materials Required but Not Supplied.....	3
Equipment	3
Supplies.....	4
Section 6. Safety Precautions and Recommendations for Best Results	5
Section 7. Protocol	7
Sample Enrichment	7
Free DNA Removal Treatment.....	8
DNA Extraction.....	8
Real-Time PCR	9
Data Analysis.....	9
Section 8. Confirmation of Positive Results	11
Section 9. Confirmation of Single Colonies Using iQ-Check Kit.....	12
Section 10. Test Performance and Validations	12
Section 11. References	13
Section 12. Revision History	13
Appendix — PCR Mix Calculation Guide	14

Section 1

Introduction

Cronobacter spp., formerly named *Enterobacter sakazakii*, is a ubiquitous pathogen that has been isolated from food, environmental, and clinical sources. It is associated with rare but fatal infant infections linked to the consumption of reconstituted powdered infant formula. Low birth weight and immunocompromised newborns, as well as neonates (infants less than four weeks old), are particularly at risk because of the high mortality rate (up to 80%) but the largest number of cases occurs in older populations.

The resistance of the organism to desiccation, the low cell numbers that have been reported to cause disease, and the high volume of baby milk product sales have led many authorities to monitor for the pathogen. For example, the EU regulation requires absence of the bacteria in 10 g of sample. To comply with regulations and to fulfill Hazard Analysis and Critical Control Points (HACCP) programs, industrial companies have stricter targets and need to select highly sensitive and specific methods. The classical microbiological methods offer standardized results. They are, however, long and tedious.

iQ-Check *Cronobacter* spp. Kit is a simple and rapid qualitative test allowing the detection of specific *Cronobacter* spp. DNA sequences in infant formula and environmental samples. Up to 94 samples can be processed to obtain a reliable result in less than 24 hr. This allows a more flexible workflow and a higher throughput to implement effective measures to improve safety and hygiene.

Section 2

iQ-Check *Cronobacter* spp. Technology

The iQ-Check *Cronobacter* spp. Kit is a test based on gene amplification and detection by real-time PCR. Ready-to-use PCR reagents contain DNA primers and a DNA probe specific for *Cronobacter* spp., as well as DNA polymerase and nucleotides. Detection and data analysis are optimized for use with a Bio-Rad real-time PCR instrument, such as the CFX96 Touch Deep Well or CFX Opus Deepwell Real-Time PCR Detection System.

PCR is a powerful technique used to generate many copies of target DNA. During the PCR reaction, several cycles of heating and cooling denature the DNA. Primers then anneal to the target region, where the DNA polymerase uses the primers and deoxynucleotide triphosphates (dNTPs) to extend the DNA, creating copies of the target DNA. These copies are called amplicons.


In real-time PCR, specific probes detect the DNA during amplification by hybridizing to the amplicons. These probes are linked to a fluorophore that fluoresces only when hybridized to the target sequence. FAM is the fluorophore linked to the probe in this kit that hybridizes to the

Section 3 Kit Components

Cronobacter spp.–specific DNA sequence. In the absence of target DNA, no fluorescence will be detected. As the number of amplicons increases with each round of amplification, fluorescence intensity also increases. At the annealing step of each PCR cycle, the optical module or detector measures this fluorescence. The software associated with the instrument plots fluorescence intensity against number of cycles. This method allows a simple determination of the presence or absence of *Cronobacter* spp. in a sample.

A synthetic DNA internal control is included in the reaction mix. This control is amplified with a specific probe at the same time as the *Cronobacter* target DNA sequence and is detected by a second fluorophore. This allows for the validation of any negative result.

The iQ-Check test can detect *Cronobacter* spp. in all infant formula, infant cereals, and environmental samples previously enriched by culture in buffered peptone water with or without supplement. See Section 7 for a list of matrices associated with all validations. The test includes five main steps:

* Please refer to the user guide of the iQ-Check Free DNA Removal Solution (#10000058391) for conditions of use.

Section 3

Kit Components

The iQ-Check *Cronobacter* spp. Kit contains sufficient reagents for 96 tests (94 samples).

Reagent ID	Reagent	Quantity Provided, ml
A	Lysis II reagent	1 bottle, 20
B	Fluorescent probes	1 tube, 0.55
C	Amplification mix	1 tube, 4.4
D	PCR negative control	1 tube, 0.5
E	PCR positive control	1 tube, 0.25

Section 4

Shelf Life and Storage

Once received, the kit must be stored at 2–8°C. Reagents stored at this temperature can be used until the expiration date indicated on the tubes.

Section 5

Materials Required but Not Supplied

Equipment

- Lab paddle blender for homogenizing test samples
- Incubator for sample microbiological enrichment
- Specific for extraction in sterile 1.5 ml conical screwcap tubes
 - Benchtop centrifuge (capable of 10,000–12,000 x g)
 - Dry heat block at $37 \pm 2^\circ\text{C}$ and/or 95–100°C
 - Cell disruptor, such as a Disruptor Genie (Scientific Industries, Inc.)
- Specific for extraction in deep well plate
 - Heating thermoshaker* capable of maintaining $37 \pm 2^\circ\text{C}$ and/or 95–100°C, with a mixing speed of at least 1,300 rpm
- Vortexer
- Magnetic stir plate
- 20, 200, and 1,000 μl micropipets
- Tips for repeat pipettors; sterile, individually packaged
- Bio-Rad real-time PCR system,* for example, the CFX96 Touch Deep Well (catalog #3600037) or CFX Opus Deepwell (catalog # 17007991) Real-Time PCR Systems
- Bio-Rad iQ-Check Prep System for automated DNA extraction and PCR plate setup (catalog #3594911)

Note: We recommend using an uninterrupted power supply (UPS) with the thermal cycler and iQ-Check Prep Systems.

* Contact Bio-Rad Technical Support for information on recommended instruments.

Supplies

- Enrichment medium: Buffered Peptone Water (for example, BPW Plus, catalog #3564684, dehydrated, 500 g; 3554179, 225 ml x 6 bottles; 3555790, 5 L x 2 bags; 3555795, 3 L x 4 bags; or BPW Standard, catalog #12013259, dehydrated, 500 g; 12013258, dehydrated, 5 kg; 12013260, 5 L x 2 bags)
- Selective supplement: Vancomycin (catalog #3564145, 5 g or equivalent)
- PIF Supplement (catalog #12013322, 2 g)
- iQ-Check Free DNA Removal Solution (catalog #3594970)
- Specific for environmental samples
 - Environmental sponges
 - Environmental swabs
 - Neutralizing broth for sponges and swabs, such as Dey-Engley (D/E) or HiCap Neutralizing Broth, or Letheen broth
- Specific for extraction in tubes
 - 1.5 ml conical screwcap tubes, sterile (catalog #17010678)
- Specific for extraction in a deep well plate
 - 96-well deep well plate (iQ-Check Deep Well Microplates, catalog #3594900)
 - Plastic sealing film (TeSeE NSP Plastic Sealing Film, catalog #3590139)
 - PCR plate sealing film (X-Pierce Films, catalog #3593977, or Pre-Pierced Plate Sealing Film, #3600040, North America only)
- Specific for iQ-Check Prep System
 - 60 ml dilution container (catalog #3594904)
 - Filter tips (catalog #12014486, 50 μ l x 5,760; #12014483, 1,000 μ l x 3,840)
 - PCR mix tubes (catalog #12016673, 5 ml x 25)
- RAPID'Sakazakii Agar (catalog #3563971, 90 mm x 20 dishes; 3564976, dehydrated, 500 g)
- PCR plates, tubes, sealing tape, and caps
- Sterile filter tips adaptable to 20, 200, and 1,000 μ l micropipets
- Tips for Combitip pipets or equivalent repeat pipettors; sterile, individually packaged
- 1 and 10 ml pipets
- 2 and 5 ml sterile test tubes

Section 6 Safety Precautions and Recommendations for Best Results

- Powder-free gloves
- Distilled sterile water
- Bleach, 5%
- Cleaning agent such as DNA AWAY or RNase AWAY

Section 6

Safety Precautions and Recommendations for Best Results

- This test must be performed by trained personnel
- Pregnant women, children, the elderly, and immunocompromised individuals should not perform this method
- Samples and enrichment cultures must be handled as potentially infectious material and discarded according to local rules and regulations
- All potentially infectious material should be autoclaved before disposal
- The quality of results depends on strict compliance with Good Laboratory Practices (for example, the EN ISO 7218 standard), especially concerning PCR:
 - Never circulate laboratory equipment (pipets, tubes, etc.) from one workstation to another
 - Always use a positive control and a negative control for each series of amplification reactions
 - Do not use reagents after their expiration date
 - Vortex reagents from the kit before using them to ensure homogeneity
 - Periodically verify the accuracy and precision of pipets, as well as correct functioning of the instruments
 - Change gloves often, especially if you suspect they are contaminated
 - Clean work spaces periodically with 5% bleach and other decontaminating agents, such as DNA AWAY
 - Use powder-free gloves and avoid fingerprints and writing on tube caps. Both will interfere with data acquisition
- It is strongly advised to follow the general requirements described in the standard EN ISO 22174:2024 (Microbiology of food and animal feeding stuffs — Polymerase chain reaction (PCR) for the detection of food pathogens — General requirements and definitions)

Section 6 Safety Precautions and Recommendations for Best Results

- iQ-Check *Cronobacter* spp. Kit
 - All substances or mixtures in the test kit are classified products, according to the Globally Harmonized System (GHS). Contact with acids may cause release of toxic gases. No special precautions are necessary if used correctly. If the product is inhaled, supply fresh air and consult a doctor in case of complaints. After eye contact with the product, rinse opened eye for several minutes under running water. If the products are swallowed, induce vomiting and call for medical help
- iQ-Check Prep System
 - Improper use of the iQ-Check Prep System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the iQ-Check Prep System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of the instrument must be performed only by Bio-Rad field service engineers
- CFX96 Touch Deep Well or CFX Opus Deepwell Real-Time PCR Detection System
 - Improper use of the CFX96 Touch or CFX Opus Deep Well Real-Time PCR Detection System may cause personal injury or damage to the instrument. Some components may pose a risk of personal injury due to excessive heat if improperly handled. For safe use, the CFX96 Touch or CFX Opus Deep Well Real-Time PCR Detection System must be operated only by qualified laboratory personnel who have been appropriately trained. Servicing of instrument must be performed only by Bio-Rad Field Service Engineers
- Enrichment
 - The user should read, understand, and follow all safety information in the instructions for the iQ-Check *Cronobacter* spp. Kit. Retain the safety instructions for future reference. To reduce the risks associated with exposure to chemicals and biohazards, perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples. Avoid contact with the contents of the enrichment media and reagent tubes after amplification. Dispose of enriched samples according to current industry standards
 - *Cronobacter* spp. are Biosafety Level 2 organisms. Biological samples such as enrichments have the potential to transmit infectious diseases. Follow all applicable local, state/provincial, and/or national regulations on disposal of biological waste. Wear appropriate protective equipment, which includes but is not limited to: protective eyewear, face shield, clothing/lab coat, and gloves. All work should be conducted in properly equipped facilities utilizing the appropriate safety equipment (for example, physical containment devices). Individuals should be trained in accordance with applicable regulatory and company/institution requirements before working with potentially infectious materials

Section 7 Protocol

- When testing is complete, all materials and media possibly containing pathogens should be decontaminated following current industry standards for the disposal of contaminated waste (that is, autoclave for 20 min at 120°C). Consult the Safety Data Sheet for additional information and local regulations for disposal

Section 7

Protocol

A. Sample Enrichment

It is strongly recommended to read the entire protocol before starting the test.

Enrichment media must be at the appropriate incubation temperature (ambient or 37°C when required) before use. The use of alpha-amylase is recommended for cereals or starch products enrichment as described in the ISO 6887-4 standard.

The following table outlines the different protocols that can be used, depending on the application and the scope of the validation.

NF VALIDATION BRD 7/23-01/13		
Scope (Matrices)	Sample Preparation	Enrichment/DNA Extraction
Environmental samples ^{1,2,3}	Homogenize <i>n</i> g sample in 9 x <i>n</i> ml BPW	Easy I extraction <ul style="list-style-type: none">▪ Incubate 16-20 hr at 37 ± 1°C▪ Tube/deep well format
Infant formula and cereal w/ and w/o probiotics (30 g) ^{1,2,3}	Homogenize <i>n</i> g food in 9 x <i>n</i> ml BPW + vancomycin then subculture in BPW	Easy I extraction <ul style="list-style-type: none">▪ Incubate BPW + vancomycin 18-22 hr at 37 ± 1°C▪ Subculture in BPW 3-5 hr at 37 ± 1°C▪ Tube/deep well format
Infant formula and cereals w/ and w/o probiotics, ingredients (up to 375 g) ^{1,4}	Homogenize <i>n</i> g food in 3 x <i>n</i> ml prewarmed BPW + PIF Supplement	Easy I extraction <ul style="list-style-type: none">▪ Incubate 18-26 hr at 37 ± 1°C▪ Tube/deep well format

¹ Validation also includes the use of the iQ-Check Free DNA Removal Solution and of the "Crono Fast" Application Protocol File for a reduced PCR run time. Contact your Bio-Rad representative for more information

² Within the scope of the NF VALIDATION, test portions weighing more than 30 g have not been tested

³ The enriched samples can be stored at 2-8°C for 48 hr following the end of the incubation at 37°C. Standard DNA extraction protocol must then be applied. Contact your Bio-Rad representative for more information

⁴ The enriched samples can be stored at 2-8°C for 72 hr following the end of the incubation at 37°C

B. Free DNA Removal Treatment

The iQ-Check Free DNA Removal Solution provides an ideal way to remove free DNA. Follow Bio-Rad's recommendations in the user guide.

C. DNA Extraction

General recommendations:

- Turn on the heat block or thermoshaker to preheat before starting the test. Set it to 95–100°C. Keep the lysis reagent in suspension while pipetting by stirring at medium speed on a magnetic stir plate
- In general, avoid shaking the enrichment bag and collecting large fragments of food debris. For food samples with a fatty supernatant, collect the sample just below this layer
- Open tubes and wells carefully to avoid any possible cross contamination
- Cool the deep well plate before pipetting directly through pre-pierced sealing film
- Use the magnetic bar to keep the lysis reagent in suspension. Pipet while it is stirring at medium speed
- Gently shake the lysis reagent by hand first to resuspend the resin. Then pipet while the magnetic stir bar is stirring the bottle contents at medium speed in order to keep the contents in suspension

Easy I Protocol

1. Aliquot 100 µl of homogenized lysis reagent (reagent A) to tubes or wells of a deep well plate.

Note: Gently shake the lysis reagent by hand first to resuspend the beads.

2. Add 100 µl of enriched sample.

Note: Shake the suspension to homogenize the culture and then allow any debris to settle before collecting the sample.

3. Mix the solution by pipetting up and down until homogenized.
4. Close the tubes or seal the deep well plate with pre-pierced sealing film.
5. Incubate tubes in the heat block at 95–100°C for 15 - 20 min. Incubate deep well plate in the thermoshaker under agitation at 1,300–1,600 rpm at 95–100°C for 15–20 min.
6. Vortex tubes at high speed, and then centrifuge at 10,000–12,000 x g for at least 2 min. Centrifugation is not needed for deep well plate.

This is the recommended stopping point for temporarily stopping the procedure.

The supernatant can be stored for up to 1 year at -20°C . Always allow it to thaw and homogenize, and then centrifuge at 10,000–12,000 x g for 5 min before reusing.

D. Real-Time PCR

Instrument and Software Setup

For instrument and software setup, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

PCR Mix Preparation

1. Prepare PCR mix containing the amplification solution (reagent C) and the fluorescent probes (reagent B). Volume of PCR mix needed depends on the number of samples and controls to be analyzed. At least one positive and one negative control must be included in each PCR run. Use the pipetting table in the Appendix to find the correct volumes to use for each reagent.
2. Use the PCR mix (reagents B + C) immediately after preparation. It is stable for 1 hr maximum at $2-8^{\circ}\text{C}$.
3. Pipet 45 μl of PCR mix into each well according to your plate setup.
4. Add 5 μl of DNA extract, reagent D (negative control), or reagent E (positive control). Do not vortex the sample before pipetting. Hermetically seal the wells of the plate or strips. It is important to avoid bubbles at the bottom of the wells by pipetting carefully. As an optional step to eliminate any bubbles, centrifuge the sealed PCR plate or the PCR strips (quick spin).
5. Place the plate or strips in the thermal cycler. Be sure to place the plate with the A1 well at the upper left corner. Close the reaction module.

Run PCR

To start the PCR run, follow instructions in the real-time PCR system user guide for iQ-Check Kits.

The Application Protocol File "Crono Fast" has been NF Validated on production environmental samples category only.

E. Data Analysis

Data can be analyzed directly at the end of the PCR run or at a later time by opening the stored data file. Follow instructions in the corresponding CFX Manager IDE Software or CFX Maestro IDE Software User Manual for opening data files and setting the data analysis parameters.

Interpreting Results

Once the data analysis parameters have been set, results are interpreted by analyzing the C_q values of each sample (the cycle at which the amplification curve crosses the threshold).

Section 7 Protocol

CFX Manager IDE Software or CFX Maestro IDE Software allows complete automated analysis for Bio-Rad real-time PCR detection systems. Verification of the typical characteristics of the amplification curves should be performed prior to releasing results. Contact your local Bio-Rad technical support team for additional support.

Controls

Verify the positive and negative controls before interpreting sample results.

For the experiment to be valid, the controls must have the following results, as summarized in the table below. Otherwise the PCR reaction must be repeated.

	<i>Cronobacter</i> spp. Detection (FAM channel)	Internal Control Detection (HEX channel)
Negative control	Cq = N/A*	$28 \leq Cq \leq 40$
Positive control	$26 \leq Cq \leq 36$	N/A

* The software indicates a Cq value (the cycle at which the amplification curve crosses the threshold) of N/A (not applicable) when the fluorescence of a sample does not rise significantly above the background noise and hence does not cross the threshold.

If results of negative and positive controls differ from those in the table above (invalid control), repeat the run and analysis described in D. Real-Time PCR and E. Data Analysis in Section 7 Protocol

Samples

A positive iQ-Check *Cronobacter* spp. PCR test must show a typical amplification curve and a Cq value ≥ 10 for the FAM fluorophore.

- If the Cq value for both channels is below 10, verify that the raw data the curve is a regular amplification curve (with a flat baseline followed by a rapid exponential increase of fluorescence and then a flattening out). If the curve seems correct, it may be considered a positive *Cronobacter* spp. test

If there is no Cq value (Cq = N/A) for FAM, or if the curve is not a typical amplification curve, the internal control for that sample must be analyzed:

- This sample is considered a negative *Cronobacter* spp. sample if there is no Cq value in the FAM channel and the internal control has a Cq value ≥ 28
- Should the internal control also not have a Cq value (Cq = N/A), this probably indicates an inhibition of the PCR reaction. The sample needs to be diluted (using 10 μ l of DNA extract, perform a 1:10 dilution in distilled sterile water, and then test 5 μ l of the dilution) and the PCR repeated
- Should the Cq value for the internal control be < 28 , it is not possible to interpret the result. Verify that the threshold was correctly placed, or that the curve as raw data is a regular

Section 8 Confirmation of Positive Results

amplification curve. If the curve does not have a characteristic shape, it will be necessary to repeat the PCR test

Interpretation of sample results is summarized in the following table:

<i>Cronobacter</i> spp. Detection (FAM channel)	Internal Control Detection (HEX channel)	Interpretation
Cq ≥ 10	N/A	Positive
Cq = N/A	Cq ≥ 28	Negative
Cq = N/A	Cq = N/A	Inhibition*

* When both target and internal control detection give a Cq value = N/A, the sample must be diluted 1:10 and tested again.

An invalid interpretation can be given when validation criteria are not met. Check the raw data and proceed as if the sample was inhibited.

Section 8

Confirmation of Positive Results

In the context of the NF VALIDATION certified method, all positive iQ-Check Kit results must be confirmed in the following ways:

- Using classic tests described in the standardized methods ISO 22964:2017 by going back to the sample
- Using RAPID'Sakazakii chromogenic medium. Refer to the confirmation method in the RAPID'Sakazakii Agar user guide (document #10000127307)
 - For environmental samples, use the RAPID'Sakazakii method. Sub-culture 0.1 ml of sample in 10 ml of mLST prior to streaking on RAPID'Sakazakii Agar
 - For infant formula and cereal samples with and without probiotics (30 g), streak 10 µl of last enriched sample enrichment onto RAPID'Sakazakii Agar
 - For infant formula and cereal samples using PIF supplement (375 g), streak a 10 µl loop of the sample enriched in BPW + PIF supplement on RAPID'Sakazakii Agar
- Using any other method certified by NF VALIDATION based on a principle different from that used in the iQ-Check *Cronobacter* spp. PCR test. The validated protocol of this second method must be followed entirely
- In the case of discrepant results between the iQ-Check *Cronobacter* spp. Kit and any of the confirmation options listed above, follow necessary steps to ensure valid results
- The enriched BPW + PIF supplement may be stored at 2–8°C for 72 hr maximum following the incubation at 37°C before carrying out the confirmation.

Section 9

Confirmation of Single Colonies Using iQ-Check Kit

iQ-Check *Cronobacter* spp. Kit may also be used for confirming single isolated *Cronobacter* spp. colonies on agar plates. This is officially validated by AFNOR Certification on RAPID[®] *Sakazakii* Agar.

1. Pick an isolated colony from a selective or nonselective agar plate with a toothpick, sterile loop, or other adapted consumable (for example, a pipet tip).
2. Resuspend the colony in 100 µl tryptone salt or distilled sterile water in a microcentrifuge tube. Homogenize using a vortexer.
3. Use 5 µl of the suspension with 45 µl of PCR mix (see D. Real-Time PCR, in Section 7 Protocol) and follow the rest of the iQ-Check *Cronobacter* spp. protocol for the data and result interpretation.

Section 10

Test Performance and Validations

The iQ-Check *Cronobacter* spp. Kit is specific for the *Cronobacter* genus.

NF Validation

BRD 07/23 – 01/13

ALTERNATIVE ANALYTICAL
METHODS FOR AGRIBUSINESS

<http://nf-validation.afnor.org>

iQ-Check *Cronobacter* spp. is certified NF VALIDATION as an alternative method to the reference method EN ISO 22964 (2017), for the detection of *Cronobacter* spp. in infant formulae and cereals and environmental samples. The validation followed the protocol of the NF EN ISO 16140: part 2 2016 standard and includes the use of the CFX96 and CFX Opus Deepwell Real-Time PCR Detection Systems. The use of the iQ-Check Free DNA Removal Solution is validated. The associated software is the CFX Manager IDE Software (V2.2 and later) or CFX Maestro IDE Software (version 4.0 and later). The "Crono Fast" APF is also validated for all samples. Certificate number: BRD 07/23 – 01/13. Valid until: refer to the certificate available on the AFNOR Certification website.

Section 11

References

ISO 7218:2024. Microbiology of the food and animal feeding stuffs – General requirements and guidance for microbiological examination

ISO 22964. Microbiology of the food chain – Horizontal method for the detection of *Cronobacter* spp.

ISO 16140-2. Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method

ISO 6887-4. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 4: Specific rules for the preparation of miscellaneous products

Section 12

Revision History

Release date	Document number	Change
May 2020	10000128585 Ver A	<ul style="list-style-type: none"> - Renewal and extension of AFNOR validation for: 375 g samples FDRS use as an option "Crono Fast" APF - New document design and update of references and content - Document number change – previous version 808473 REV D
December 2022	10000128585 Ver B	<ul style="list-style-type: none"> - AFNOR extension: CFX Opus Deepwell Real-Time PCR System and upgrade to CFX Manager Software IDE version 3.1 - Clarification of result interpretation protocol - Added new catalog numbers for consumables - Replaced protocol section with summary table - General content updates
December 2023	10000128585 Ver C	<ul style="list-style-type: none"> - Updated protocol table and confirmation method
December 2024	10000128585 Ver D	<ul style="list-style-type: none"> - Clarification of the infant formula and cereal w/ and w/o probiotics protocol
April 2025	10000128585 Ver E	<ul style="list-style-type: none"> - Validation extension for CFX Maestro Software IDE version 4.0 - Added the option to use an equivalent reference for vancomycin

Appendix — PCR Mix Calculation Guide

To find the correct volumes to use when preparing the PCR mix, add the total number of samples and controls to be analyzed and find the corresponding volumes of reagent B and reagent C in the table.

Total Number of Samples and Controls	Probes Reagent B, μl	Amplification Mix Reagent C, μl	Total Number of Samples and Controls	Probes Reagent B, μl	Amplification Mix Reagent C, μl	Total Number of Samples and Controls	Probes Reagent B, μl	Amplification Mix Reagent C, μl
1	5	40	33	178	1,400	65	351	2,800
2	11	86	34	184	1,500	66	356	2,900
3	16	130	35	189	1,500	67	362	2,900
4	22	173	36	194	1,600	68	367	2,900
5	27	216	37	200	1,600	69	373	3,000
6	32	259	38	205	1,600	70	378	3,000
7	38	302	39	211	1,700	71	383	3,100
8	43	346	40	216	1,700	72	389	3,100
9	49	389	41	221	1,800	73	394	3,200
10	54	432	42	227	1,800	74	400	3,200
11	59	475	43	232	1,900	75	405	3,200
12	65	518	44	238	1,900	76	410	3,300
13	70	562	45	243	1,900	77	416	3,300
14	76	605	46	248	2,000	78	421	3,400
15	81	648	47	254	2,000	79	427	3,400
16	86	691	48	259	2,100	80	432	3,500
17	92	734	49	265	2,100	81	437	3,500
18	97	778	50	270	2,200	82	443	3,500
19	103	821	51	275	2,200	83	448	3,600
20	108	864	52	281	2,200	84	454	3,600
21	113	907	53	286	2,300	85	459	3,700
22	119	950	54	292	2,300	86	464	3,700
23	124	994	55	297	2,400	87	470	3,800
24	130	1,000	56	302	2,400	88	475	3,800
25	135	1,100	57	308	2,500	89	481	3,800
26	140	1,100	58	313	2,500	90	486	3,900
27	146	1,200	59	319	2,500	91	491	3,900
28	151	1,200	60	324	2,600	92	497	4,000
29	157	1,300	61	329	2,600	93	502	4,000
30	162	1,300	62	335	2,700	94	508	4,100
31	167	1,300	63	340	2,700	95	513	4,100
32	173	1,400	64	346	2,800	96	518	4,100

Visit bio-rad.com/iqcheck for more information.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

IQ-CHECK is a trademark of Bio-Rad Europe GMBH in certain jurisdictions.

All trademarks used herein are the property of their respective owner.

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ISO 22964 Microbiology of the food chain – Horizontal method for the detection of *Cronobacter spp.*

ISO 16140-2: Microbiology of the food chain – Method validation – Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method

ISO 6887-4 Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 4: Specific rules for the preparation of miscellaneous products

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