Listeria Special Broth II (LSB II)

Catalog #	Description
12020718	LSB II , 5 L x 2 bags
12017463	LSB II , 225 ml x 6 bottles
12017388	LSB II , dehydrated, 500 g
12017378	LSB II , dehydrated, 5 kg

For laboratory use only.

Intended Use

Selective broth used for the enrichment of *Listeria monocytogenes* and other *Listeria* species in food products and environmental samples when using iQ-Check Real-Time PCR Detection Kits.

Principle

The peptone mixture, yeast extract and nutritive mixture provide nutritional elements necessary for the growth of *Listeria*. The presence of an adequate buffer ensures control of pH and optimizes recovery of injured organisms. The selective mixture inhibits competing microorganisms.

Formula

Peptone mixture	15 g
Yeast extract	5 g
Buffer	22.2 g
Ferric citrate (Fe III)	0.5 g
Selective mixture	8 g
Nutritive mixture	3 g
Distilled water	1,000 ml
Final pH at 25°C =	7.2 ± 0.2

Shelf Life and Storage

Store dehydrated media at 15–25°C in a carefully sealed package in a dry and dark place. Store the prepared liquid media in sterile conditions and a dark place at 2–8°C for 1 month or at 12–25°C for 1 week. Store liquid media at 2–8°C until the expiration date or up to 3 months at room temperature (30°C maximum) in a dark place.

Required Materials Not Supplied

This is a non-exhaustive list.

Equipment

- All usual laboratory equipment
- Incubators or incubation room
- Scales
- Stirrer/homogenizer
- Vortexer

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- The medium may produce a white residue on the walls of the bottle. This has no effect on quality
- Bags of broth must be manipulated under aseptic conditions. After opening, the bags can be stored carefully sealed with a sterile clamp at ambient temperature for 3 days. Bags can be re-opened once. Do not use the volume remaining in the spout, remove it before dispensing the broth
- If the necessary, LSB II broth can be pre-heated at 37°C before use. For laboratory organization, the broth could be maintained at 37°C up to 1 week
- Some of the components of the dehydrated product represent a hazard. Refer to the Safety Data Sheet (SDS)
- For SDS product safety information and certificate of analysis, visit bio-rad.com

Quality Control

Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.

Protocol

Dehydrated LSB II Preparation

- Always shake the bottle before use
- Dissolve 53.7 g of LSB II in 1 L of sterile distilled water
- Mix until a homogeneous suspension is obtained
- Dispense in appropriate container
- Sterilize in autoclave at $110 \pm 3^\circ\text{C}$ for 15 min

Sample Preparation and Enrichment Protocol

- Prepare sample according to the standard method applicable to the product concerned
- Refer to the user guides for the iQ-Check *Listeria monocytogenes* (document #10000167776) and iQ-Check *Listeria* spp. (document #10000167777) Kits for enrichment protocols

References

- Curtis GDW et al. (1989). A selective differential medium for the isolation of *Listeria monocytogenes*. Lett Appl Microbiol 8, 95–98.
- Curtis GDW et al. (1989). Selective agents for *Listeria* can inhibit their growth. Lett Appl Microbiol 8, 169–172.
- ISO 7218:2007/AMD1:2013. Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.
- Lovett J et al. (1987). *Listeria monocytogenes* in raw milk: Detection, incidence, and pathogenicity. J Food Prot 50, 188–192.
- Prentice GA and Neaves P (1988). *Listeria monocytogenes* in Food: Its Significance and Methods for Its Detection. (Brussels: International Dairy Federation), p. 223.
- Van Netten P et al. (1988). A selective and diagnostic medium for use in the enumeration of *Listeria* spp. in foods. Int J Food Microbiol 6, 187–198.

Revision History

Release date	Document number	Change
January 2024	10000171166 Ver A	<ul style="list-style-type: none">- New document number – previous version 10000154288 Ver B <i>Listeria</i> Special Broth II User Guide- Addition of liquid bags format- Storage conditions update- Precautions update

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. IQ-CHECK is a trademark of Bio-Rad Europe GMBH in certain jurisdictions. All trademarks used herein are the property of their respective owner.

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Riferimenti

Curtis GDW et al. (1989). A selective differential medium for the isolation of *Listeria monocytogenes*. Lett Appl Microbiol 8, 95–98.

Curtis GDW et al. (1989). Selective agents for *Listeria* can inhibit their growth. Lett Appl Microbiol 8, 169–172.

ISO 7218:2007/AMD1:2013. Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.

Lovett J et al. (1987). *Listeria monocytogenes* in raw milk: Detection, incidence, and pathogenicity. J Food Prot 50, 188–192.

Prentice GA and Neaves P (1988). *Listeria monocytogenes* in Food: Its Significance and Methods for Its Detection. (Brussels: International Dairy Federation), p. 223.

Van Netten P et al. (1988). A selective and diagnostic medium for use in the enumeration of *Listeria* spp. in foods. Int J Food Microbiol 6, 187–198.

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Curtis GDW et al. (1989). A selective differential medium for the isolation of *Listeria monocytogenes*. Lett Appl Microbiol 8, 95 – 98.

Curtis GDW et al. (1989). Selective agents for *Listeria* can inhibit their growth. Lett Appl Microbiol 8, 169 – 172.

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Lovett J et al. (1987). *Listeria monocytogenes* in raw milk: Detection, incidence, and pathogenicity. J Food Prot 50, 188 – 192.

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Van Netten P et al. (1988). A selective and diagnostic medium for use in the enumeration of *Listeria* spp. in foods. Int J Food Microbiol 6, 187 – 198.

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