
RAPID'*P.aeruginosa* Agar

User Guide

Chromogenic media for the direct enumeration without confirmation of *Pseudomonas aeruginosa* in water samples using the membrane filtration method

Catalog #3563984, Prepared plates, 55 mm x 20 dishes
Catalog #3564900, Dehydrated, 500 g

BIO-RAD

Table of Contents

Section 1	Introduction	1
Section 2	RAPID' <i>P.aeruginosa</i> Principle	1
Section 3	Theoretical Formula	1
Section 4	Shelf Life and Storage	1
Section 5	Materials Required but Not Supplied	1
	Equipment.....	1
	Supplies	2
Section 6	Precautions, Limitations of Use, and Quality Control	2
Section 7	Protocol	3
	Preparation of Dehydrated Medium.....	3
	Sample Preparation	3
	Inoculation and Plate Reading.....	3
Section 8	Confirmation of Positive Results.....	3
Section 9	Confirmation of Other Methods.....	4
Section 10	Test Performance and Validations.....	4
Section 11	References.....	4
Section 12	Revision History	5

Section 1 Introduction

Pseudomonas aeruginosa is a ubiquitous environmental gram-negative bacterium naturally present in fresh and seawater, moist soil, and on the surface of vegetables. *P. aeruginosa* is an opportunistic pathogen that can cause serious illness. Many companies impose regulations on *P. aeruginosa* testing on water for human consumption, with absence as the main criteria. These regulations are often adopted where government regulations do not mandate this testing in order to ensure the safety of drinking water.

Section 2 RAPID'*P.aeruginosa* Principle

The principle of RAPID'*P.aeruginosa* medium is based on the detection of an enzymatic activity typical of *P. aeruginosa*. Under its action, a specific chromogenic substrate is cleaved, leading to the formation of a blue to blue-green or green precipitate on *P. aeruginosa* colonies. The selective mixture makes it possible to inhibit the majority of interfering flora, in particular other *Pseudomonas non-aeruginosa* strains. Other microorganisms may, however, show growth. Their colonies appear transparent or pigmented yellow-green and are easily distinguishable from those of *P. aeruginosa*.

Section 3 Theoretical Formula

Nutritive mix	14 g
Buffer system	2.15 g
Selective agents	0.11 g
Chromogenic mix	0.13 g
Agar	10 g
Distilled water	qsp 1,000 ml

Final pH at 25°C = 7.1 ± 0.2

Section 4 Shelf Life and Storage

- Dehydrated medium: 2–8°C in carefully sealed package, in a dry and dark place
- Prepoured plates: 2–8°C in a dark place
- Plate prepared from dehydrated medium: 15 days at 2–8°C in carefully sealed package, in a dry and dark place

Section 5 Materials Required but Not Supplied

Equipment

- All usual laboratory equipment
- Filtration device

- Scale, sensitivity of 0.1 g
- Stirrer/homogenizer
- Thermostatically controlled incubator or incubation room, precise to $\pm 1^{\circ}\text{C}$
- Water bath, precise to $\pm 1^{\circ}\text{C}$

Supplies

- Forceps for handling membranes
- Sterile distilled water
- Sterile filtration membranes ($\varnothing = 47$ mm, 0.45 μm Millipore HAWG 047 Type HA)
- Sterile petri dishes ($\varnothing 55$ mm)
- Sterile pipets

Section 6

Precautions, Limitations of Use, and Quality Control

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria such as *Pseudomonas*
- Media that have come in contact with water samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Avoid trapping air bubbles underneath the membrane during its placement on the agar. Poor membrane-agar contact may lead to an erroneous result. If necessary, gently and carefully flatten the membrane with the forceps

Limitations of Use

Not applicable

Quality Control

- Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file
- For SDS product safety information and certificate of analysis, visit www.bio-rad.com

Section 7 Protocol

Preparation of Dehydrated Medium

1. Always shake bottle before use.
2. Dissolve 26.4 g of powder in 1 L of distilled water and mix until a homogenous suspension is obtained.
3. Heat gently, agitating frequently, then bring to a boil until completely dissolved.
4. Sterilize in an autoclave at $121 \pm 3^\circ\text{C}$ for 15 min.
5. Cool the medium to $44\text{--}47^\circ\text{C}$. Dispense in petri dishes and leave on the bench to dry.
6. One 500 g bottle of powder makes 18.9 L of medium.

Sample Preparation

To be carried out in accordance with the standard of the product concerned.

Inoculation and Plate Reading

1. Bring medium to room temperature.
2. Under sterile conditions, filter the volume of water to be analyzed through a membrane according to the origin of the sample (e.g., 250 ml for bottled water samples or 100 ml for pool water samples).
3. Place the membrane, cross-hatched surface up, on the surface of the medium, taking care that the membrane-agar contact is complete.
4. Incubate plate at $36 \pm 2^\circ\text{C}$ for 22–30 hr.
5. *P. aeruginosa* form blue, blue-green or green colonies on RAPID'*P.aeruginosa* agar. All colonies demonstrating typical color should be counted, regardless of size.
6. Refer to ISO 8199 standard for inoculation, enumeration and interpretation of results.
7. Record the total number of *P. aeruginosa* per unit volume of filtered sample.

Section 8 Confirmation of Positive Results

Not applicable

Section 9 Confirmation of Other Methods

Not applicable

Section 10 Test Performance and Validations

Certification	Scope	Validation Protocol	Reference Protocol	Certificate Reference
NF VALIDATION	<p>Bottled water (250 ml)</p> <p>Water for human consumption (wells, springs and boreholes, tap water and drinking fountain water) (100 ml)</p> <p>Treated recreational water (pool water, thermal water) (100 ml)</p>	NF 148	EN ISO 16266	 <p>BRD 07/21-04/12 Water analysis methods http://nf-validation.afnor.org/en</p> <p>For more information about the end of validity of the NF VALIDATION certification, please refer to the certificate BRD 07/21-04/12 available on the website http://nf-validation.afnor.org/en and/or on request asking Bio-Rad Laboratories</p>

Section 11 References

Balasubramanian D, Mathee K. (2009). Comparative transcriptome analyses of *Pseudomonas aeruginosa*. *Genomics* 4, 349-361.

EN ISO 11133: Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.

ISO 17025: Acceptance sampling plans and procedures for the inspection of bulk materials.

Lisa A.T., P.R. Beassoni, M.J. Massimelli, L.H. Otero, C.E. Domenech (2007). *Appl Microbiol.* 255-262.

Lucchesi G.I., C. Pallotti, A. T. Lisa, C. E. Domenech (1998). *Microbiol Letters* (12), 123-126.

EN ISO 16266: Water quality – Detection and enumeration of *Pseudomonas aeruginosa* – Membrane filtration method.

Section 12 Revision History

Release date	Document number	Change
January 2024	10000171245 Ver A	<ul style="list-style-type: none"> - Document number change – previous version 100000126751 Ver C RAPID' <i>P.aeruginosa</i> User Guide - General content updates

Visit www.bio-rad.com/water for more information on our water testing products.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

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