AL Agar

User Guide

Chromogenic media for the detection and enumeration of the *Listeria monocytogenes* and other species of *Listeria* in food products for human and animal consumption, and in environmental samples

Catalog #3563695, Prepared plates, 90 mm x 20 dishes
Catalog #3563965, Prepared plates, 90 mm x 120 dishes
Catalog #3555200, Bottled agar base, 237.5 ml x 6 bottles
Catalog #3564043, Dehydrated, 500 g
Catalog #3564041, Supplement 1, freeze dried, 10 vials
Catalog #3564042, Supplement 2, liquid 25 ml, 10 vials


Table of Contents

Section 1	Introduction	1
Section 2	AL Agar Principle	1
Section 3	Theoretical Formula.....	1
Section 4	Shelf Life and Storage	2
Section 5	Materials Required but Not Supplied	2
	Equipment.....	2
	Supplies	2
Section 6	Precautions, Limitations of Use, and Quality Control	3
Section 7	Protocol.....	4
	Preparation of Dehydrated Medium.....	4
	Preparation of Bottled Medium	4
	Detection of <i>L. monocytogenes</i> and <i>Listeria</i> Genus.....	4
	Enumeration of <i>L. monocytogenes</i>	5
Section 8	Confirmation of Positive Results	6
Section 9	Confirmation of Other Methods.....	6
Section 10	Test Performance and Validation	7
Section 11	References.....	7
Section 12	Revision History.....	8

Section 1 Introduction

Serious outbreaks of *Listeria monocytogenes* continue to plague the food safety industry. *Listeria* is dangerous in the way that it can survive and grow, however slowly, at refrigerated temperatures in ready-to-eat processed foods. *Listeria monocytogenes* is a particularly problematic pathogen, as it can cause serious health problems and possible death, particularly in immunosuppressed individuals, newborns, and the elderly. Infection is known to cause still births and miscarriages in pregnant women. Listeriosis develops with symptoms such as fever, fatigue, nausea, vomiting, and diarrhea and has a mortality rate of 30%, which may be higher in vulnerable individuals. Approximately 90% of all reported cases of listeriosis result in hospitalization. A rapid culture method is necessary to cut down on time to results and to make sure those results are accurate.

Section 2 AL Agar Principle

The principle of AL medium (Agar *Listeria* according to Ottaviani and Agosti) is based on the simultaneous detection of two enzyme activities: β -D-glucosidase and phosphatidylinositol-specific phospholipase C (PIPLC). The β -D-glucosidase activity, common to all *Listeria*, is detected using a chromogenic substrate (X-glucoside). Hydrolysis of this substrate by the *Listeria* genus leads to the production of blue to blue-green colored colonies. PIPLC is an enzyme detected only in pathogenic *Listeria* species, *L. monocytogenes* and *L. ivanovii*. AL medium contains phosphatidylinositol, which when broken down, produces an opaque halo around colonies of bacteria of these two species. This halo generally appears after 24 hr of incubation in *L. monocytogenes* and after 48 hr of incubation in *L. ivanovii*. Selectivity of the medium is achieved by the combined action of lithium chloride, antibiotics, and an antifungal. This medium may be used in standardized protocol or in an alternative validated short protocol.

Section 3 Theoretical Formula

Meat peptone	18 g
Tryptone	6 g
Yeast extract	10 g
Sodium pyruvate	2 g
Glucose	2 g
Anhydrous magnesium glycerophosphate	1 g
Anhydrous magnesium sulphate	0.5 g
Sodium chloride (NaCl)	5 g
Lithium chloride (LiCl)	10 g
Anhydrous Na ₂ HPO ₄	2.5 g
Chromogenic substrate (5-bromo-4-chloro-3-indolyl- β -D-glucopyranoside)	0.05 g
Nalidixic acid	0.02 g
Ceftazidime	0.02 g
Polymyxin B sulphate	76,700 UI
Amphotericin B	0.01 g
Phosphatidylinositol	2 g
Agar	12 g
Distilled water	qsp 1,000 ml

Final pH at 25°C = 7.2 ± 0.2

Section 4 Shelf Life and Storage

- Dehydrated: 15–25°C in carefully sealed package, in a dry and dark place
- Bottled agar base and supplements: 2–8°C in a dark place
- Reconstituted supplement 1: 1 week at 2–8°C
- Pre-poured: 2–8°C in a dark place
- Plate prepared from the dehydrated base: 2 weeks at 2–8°C in carefully sealed package, in a dry and dark place

Section 5 Materials Required but Not Supplied

Equipment

- All usual laboratory equipment
- Hot plate
- Scale, sensitivity of 0.1 g
- Stirrer/homogenizer
- Thermostatically-controlled incubator or incubation room, precise to $\pm 1^\circ\text{C}$
- Water bath

Supplies

- Confirmation:
 - iQ-Check *Listeria monocytogenes* II Real-Time PCR Detection Kit (catalog #3578124), iQ-Check *Listeria* spp. Real-Time PCR Detection Kit (catalog #3578113)
 - RAPID'*L.mono* Agar (catalog #3563694, prepared plates, 90 mm x 20 dishes; 3563964, prepared plates, 90 mm x 120 dishes; 3555294, ready-to-use kit (for 200 ml) including bottled agar medium and supplements; 3564293, dehydrated, 500 g; 3564294, supplement 1, freeze-dried, 10 vials; 3564746, supplement 2, liquid, 10 vials)
 - PALCAM agar (for example catalog #3564754, dehydrated 500 g; 3564752, supplement, 10 vials)
 - Rhamnose Test (catalog #3553669, 1 ml x 28 vials)
- Diluent for enumeration:
 - Tryptone Salt (catalog #3555754, 9 ml x 25 tubes; 3555756, 900 ml x 6 bottles; 3555796, 3 L x 4 bags; 3564544, 500 g)
 - Buffered Peptone Water (for example, BPW Plus, catalog #3564684, dehydrated, 500 g; 3554179, 225 ml x 6 bottles; 3555790, 5 L x 2 bags; 3555795, 3 L x 4 bags;

Section 6 Precautions, Limitations of Use, and Quality Control

or BPW Standard, catalog #12013259, dehydrated, 500 g; 12013258, dehydrated, 5 kg; 12013260, 5 L x 2 bags)

- Enrichment medium: Half Fraser Broth (catalog #3555797, 225 ml bottles x 6; 3555794, 3 L x 4 bags; 3564604, dehydrated 500 g; 3564616, supplement, 10 vials) or Half Fraser Standard Broth (catalog #12018707 5 L x 2 bags; 12018706 dehydrated 500 g; 12018705 dehydrated 5 kg)
- Inoculating loops and spreaders
- Sterile petri dishes (90 mm)
- Sterile pipets
- Sterile weigh bags

Section 6 Precautions, Limitations of Use, and Quality Control

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria such as *Listeria*
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Allow plates to dry at 25–50°C according to ISO standard 7218 until droplets disappear from the surface of the medium. However, avoid prolonged drying, so as to not modify the efficiency of the medium. Do not dry plates in a laminar flow hood
- For heavily-loaded dishes with intensely opaque agar, reading can be facilitated by comparing opacity of the agar with a non-inoculated AL dish
- This medium may be used in standard protocols (for example, ISO or FDA BAM)
- For optimal performance of the alternative method, the use of Bio-Rad's Half Fraser broth is recommended.

Limitations of Use

- Other gram positive β -D-glucosidase positive bacteria exist without halos (for example, *Enterococcus* spp.) and with halos (for example, *Bacillus circulans*)
- Rare strains of *Enterococcus* can, within 48 hr of reading, exhibit slow phosphatidylinositol activity and ferment rhamnose like *Listeria monocytogenes*. A catalase test, directly applied to an isolated blue colony with a halo, allows differentiation between *Listeria monocytogenes* and *Listeria ivanovii* (catalase positive) from *Enterococcus* spp. (catalase negative). More generally, the atypical characteristics of *Listeria* strains are cited in the standard NF EN ISO 11290-1 (Annex D)
- Sample volumes over 25 g have not been tested as part of NF VALIDATION

Quality Control

- Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.
- For SDS product safety information and certificate of analysis, visit www.bio-rad.com.

Section 7 Protocol

Preparation of Dehydrated Medium

1. Always shake bottle before use.
2. Dissolve 34.55 g of powder in 470 ml of distilled water and mix until a homogenous suspension is obtained.
3. Heat gently, agitating frequently, then bring to a boil.
4. Sterilize by autoclaving at 121°C for 15 min.
5. Cool the medium to 44–47°C.
6. Aseptically rehydrate lyophilized supplement 1 with 5 ml sterile distilled water. Stir until completely dissolved and add to molten agar.
7. Aseptically add 1 vial of supplement 2. The temperature of the base should not exceed 50°C when adding supplement 2.
8. Mix well and pour into petri dishes. Stack dishes together (up to 5) to allow slow cooling.
9. One 500 g bottle of powder makes 7.2 L of medium.

Preparation of Bottled Medium

1. In a boiling water bath, melt the contents of the AL Agar bottle (237.5 ml).
2. Cool the bottle to 44–47°C in a water bath.
3. Aseptically rehydrate lyophilized supplement 1 with 5 ml sterile distilled water. Stir until completely dissolved.
4. Aseptically add 2.5 ml of rehydrated supplement 1 and 12.5 ml of supplement 2.
5. Mix thoroughly. Avoid frothing.
6. Pour the medium into petri dishes. Stack dishes together (up to 5) to allow slow cooling.

Detection of *L. monocytogenes* and *Listeria* Genus

1. Dilute n g or n ml of sample in $9 \times n$ ml Half Fraser Broth. Homogenize with stirrer/homogenizer.
2. Ensure the environmental swab is completely immersed in the broth (at least 9 ml).
3. Ensure the environmental sponge, cloth or gauze pad is completely soaked in broth. Add broth at 9 times the weight of the moistened wiping device.
4. Incubate at $30 \pm 1^\circ\text{C}$ for 18–28 hr.
5. After the enrichment step, the broth can be stored at $2\text{--}8^\circ\text{C}$ for 72 hr.
6. Using a sterile pipet remove 0.1 ml of sample and place drops on the outside edge of half of the agar surface.
7. Using a sterile Pasteur pipet or inoculating loop, spread sample over half the agar surface in a to-and-fro motion.
8. On the other half of the agar surface, streak for isolation by spreading the deposit in relatively close streaks over the entire dish from the edge of the previous spread.
9. Incubate at $37 \pm 1^\circ\text{C}$ for 24 ± 2 hr. Incubation can be extended to 48 ± 2 hr.
10. *L. monocytogenes* form blue/blue-green colonies with a halo.
11. Other *Listeria* spp. form blue/blue-green colonies with or without a halo.
12. After the incubation step, plates can be stored at $2\text{--}8^\circ\text{C}$ for 72 hr.
13. In case of characteristic *L. monocytogenes* colonies when using the detection protocol, it is not necessary to perform a confirmation if the sample has been already confirmed as positive in the enumeration protocol.

Enumeration of *L. monocytogenes*

1. Dilute n g or n ml of sample in $9 \times n$ ml of Half Fraser Broth or BPW.
2. Spread 0.1 ml of sample or its decimal dilutions per plate or transfer 1 ml into an empty petri dish and pour 15 ml of melted ($44\text{--}47^\circ\text{C}$) AL agar
3. If estimating small numbers, spread 1 ml of sample or decimal dilutions over three plates of 90 mm diameter (~ 0.33 ml/plate) or over one plate of 140 mm diameter.
4. Incubate the upturned dishes at $37 \pm 1^\circ\text{C}$ for 48 ± 3 hr. A first reading at 24h allows a more rapid detection of heavily contaminated samples. However the final result count is reached after 48 h (± 3 h)
5. Read plate and enumerate typical colonies.
6. *L. monocytogenes* form blue/blue-green colonies with a halo.
7. Other *Listeria* spp. form blue/blue-green colonies with or without a halo.

Section 8 Confirmation of Positive Results

8. After the incubation step, plates can be stored at 2–8°C for 72 hr.
9. Refer to EN ISO 7218 standard for inoculation, colony counting, calculation and expression of results.
10. In case of characteristic *L. monocytogenes* colonies when using the enumeration protocol, it is not necessary to perform a confirmation if the sample has been already confirmed as positive in the detection protocol.

Section 8 Confirmation of Positive Results

1. In the context of NF VALIDATION, all positive results must be confirmed. When following the enumeration protocol, if the first colony confirmed gives a negative result, confirmations should continue on up to five colonies. Confirming fewer than five colonies involves a risk of making an overestimation because of the presence of typical colonies that are not *L. monocytogenes*.
2. Confirmation of *L. monocytogenes* must be done in one of the following ways:
 - a. Using the conventional tests described in the ISO standard methods (including the purification step).
 - b. Using nucleic probes as described in the ISO 7218 standard (for example, iQ-Check *L. monocytogenes* II Real-Time PCR Detection Kit, catalog #3578124) using isolated colonies (with or without purification step).
 - c. Confirmed using a Rhamnose Test (catalog #3553669).
 - d. Repicking and spotting at least one isolated colony on a RAPID'*L.mono* agar plate. Up to 12 colonies can be confirmed on a single RAPID'*L.mono* plate.
 - e. Using any other NF VALIDATION–certified method based on a different principle from that of AL Agar. The validated protocol of the second method must be respected in its entirety. All steps preceding the detection step used as a starting point for confirmation must be common to both methods.
3. Confirmation of *Listeria* spp. other than *L. monocytogenes* must be done in one of the following ways:
 - a. Using the conventional tests described in the ISO standard methods (including the purification step).
 - b. Using nucleic probes as described in the ISO 7218 standard (for example, iQ-Check *Listeria* spp. Real-Time PCR Detection Kit, catalog #3578113) using isolated colonies (with or without purification step).
 - c. Repicking and streaking at least one isolated colony on a PALCAM agar plate. Up to 6 colonies can be confirmed on a single PALCAM agar plate.
 - d. Using any other NF VALIDATION–certified method based on a different principle from that of AL Agar. The validated protocol of the second method must be respected in its entirety. All steps preceding the detection step used as a starting point for confirmation must be common to both methods.

Section 9 Confirmation of Other Methods

4. In the event of discordant results (presumptive positive with AL Agar, negative with confirmation method), the laboratory must follow the necessary steps to ensure validity of the result obtained.

Section 9 Confirmation of Other Methods

Not applicable.

Section 10 Test Performance and Validations

Certification	Scope	Validation Protocol	Reference Protocol	Certificate Reference
NF VALIDATION	All human food products and production environmental samples	EN ISO 16140-2	EN ISO 11290-1 EN ISO 11290-2	 BRD 07/16-01/09 BRD 07/17-01/09 ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS Certified by AFNOR Certification Valid until: Refer to the certificate available on the AFNOR Certification website: http://nf- validation.afnor.org/en

Section 11 References

ISO 7218:2024. Microbiology of the food chain - General requirements and guidance for microbiological examinations

ISO 6887-1:2017/ADM: 2024. Microbiology of the food chain — Preparation of test samples, initial suspension and decimal dilutions for microbiological examination — Part 1: general rules for the preparation of the initial suspension and decimal dilutions — Amendment 1: Requirements and guidance on the use of larger test portion size for qualitative method

ISO 18593:2018. Microbiology of the food chain - Horizontal methods for surface sampling

ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 1: Detection method.

ISO 11290-2:2017. Microbiology of the food chain — Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 2: Enumeration method.

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United States Department of Agriculture (2019). Microbiology Laboratory Guidebook. Chapter 8.11: Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry, Ready-to-Eat Siluriformes (Fish) and Egg Products, and Environmental Sponges. <https://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/mlg-8.pdf?MOD=AJPERES>, accessed August 10, 2020.

United States Food and Drug Administration. Bacteriological Analytical Manual. Chapter 10: *Listeria monocytogenes*. <https://www.fda.gov/food/laboratory-methods-food/bacteriological-analytical-manual-bam>, accessed August 10, 2020.

Vlaemynck G et al. (2000): Improvement of the detection of *Listeria monocytogenes* by the application ALOA, a diagnostic, chromogenic isolation medium. J. Appl Microbiol 88: 430-441.

Section 12 Revision History

Release date	Document number	Change
December 2024	10000245369 Ver A	<ul style="list-style-type: none">- Document number change (previous version 10000135677 Ver B)- Update of Limitations of use for rare <i>Enterococcus</i> strains- Update Reference section

Visit www.bio-rad.com/rapidmedia for more information on our complete range of RAPID chromogenic media.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc.

IQ-CHECK is a trademark of Bio-Rad Europe GMBH in certain jurisdictions.

All trademarks used herein are the property of their respective owner.

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ISO 6887-1:2017/ADM: 2024. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: general rules for the preparation of the initial suspension and decimal dilutions - Amendment 1: Requirements and guidance on the use of larger test portion size for qualitative method

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United States Department of Agriculture (2019). Microbiology Laboratory Guidebook. Chapter 8.11: Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry, Ready-to-Eat Siluriformes (Fish) and Egg Products, and Environmental Sponges.

<https://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/mlg-8.pdf?MOD=AJPERES>, accessed August 10, 2020.

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Vlaemyneck G et al. (2000): Improvement of the detection of *Listeria monocytogenes* by the application ALOA, a diagnostic, chromogenic isolation medium. J. Appl Microbiol 88: 430-441.

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ISO 6887-1:2017/ADM: 2024. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1 : general rules for the preparation of the initial suspension and decimal dilutions - Amendment 1: Requirements and guidance on the use of larger test portion size for qualitative method

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ISO 11290-1:2017. Microbiology of the food chain — Horizontal method for the detection and

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enumeration of *Listeria monocytogenes* and of *Listeria* spp. — Part 1: Detection method.

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Vlaemyneck G et al. (2000): Improvement of the detection of *Listeria monocytogenes* by the application ALOA, a diagnostic, chromogenic isolation medium. *J. Appl Microbiol* 88: 430-441.

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