---

# RAPID'*L.mono* Agar

## User Guide

**Selective chromogenic medium for detection and enumeration of *Listeria monocytogenes* and other species of *Listeria* in food products for human consumption and in environmental samples**

Catalog #3563694, Prepared plates, 90 mm x 20 dishes


Catalog #3563964, Prepared plates, 90 mm x 120 dishes

Catalog #3555294, Ready-to-use kit (for 200 ml), includes bottled agar medium, supplements

Catalog #3564293, Dehydrated, 500 g

Catalog #3564294, Supplement 1, freeze dried, 10 vials

Catalog #3564746, Supplement 2, liquid, 10 vials

**BIO-RAD**

# Table of Contents

Section 1	Introduction .....	1
Section 2	RAPID' <i>L.mono</i> Principle.....	1
Section 3	Theoretical Formula.....	1
Section 4	Shelf Life and Storage.....	2
Section 5	Materials Required but Not Supplied .....	2
	Equipment .....	2
	Supplies .....	2
Section 6	Precautions, Limitations of Use, and Quality Control .....	3
Section 7	Protocol .....	4
	Preparation of Dehydrated Medium.....	4
	Preparation of Bottled Medium.....	4
	Detection of <i>L. monocytogenes</i> and <i>Listeria</i> genus.....	5
	Enumeration of <i>L. monocytogenes</i> .....	6
Section 8	Confirmation of Positive Results .....	6
Section 9	Confirmation of Other Methods.....	7
Section 10	Test Performance and Validations.....	8
Section 11	References .....	8
Section 12	Revision History .....	10

## Section 1 Introduction

Serious outbreaks of *Listeria monocytogenes* continue to plague the food safety industry. *Listeria* is dangerous in the way that it can survive and grow, however slowly, at refrigerated temperatures in ready-to-eat processed foods. *Listeria monocytogenes* is a particularly problematic pathogen, as it can cause serious health problems and possible death, particularly in immunosuppressed individuals, newborns, and the elderly. Infection is known to cause still births and miscarriages in pregnant women. Listeriosis develops with symptoms such as fever, fatigue, nausea, vomiting, and diarrhea and has a mortality rate of 30%, but may be higher in vulnerable individuals. Approximately 90% of all reported cases of listeriosis result in hospitalization. A rapid culture method is necessary to cut down on time to results and to make sure those results are accurate.

## Section 2 RAPID'L.mono Principle

The principle behind RAPID'L.mono Medium relies on the specific detection of phospholipase C (PIPLC) activity of *L. monocytogenes* and on the inability of this species to metabolise xylose. After 24 hr of incubation, *L. monocytogenes* form characteristic blue (pale blue, grey-blue to dark blue) colonies without a yellow halo. Colonies formed by other species of *Listeria* are white, with or without a yellow halo. *Listeria ivanovii* form blue-green colonies with a yellow halo (xylose-positive characteristic). This halo can appear after 24–48 hr of incubation. The selective solution in the medium permits inhibition of most interfering flora (gram-positive and gram-negative bacteria, yeasts, and molds). RAPID'L.mono permits a rapid and specific identification of *L. monocytogenes* in 24 hr and of other *Listeria* species in 24–48 hr after preparing samples in compliance with standards.

## Section 3 Theoretical Formula

Peptones	30 g
Meat extract	5 g
Yeast extract	1 g
Lithium chloride	9 g
Xylose	10 g
Phenol red	0.12 g
Growth activators	2 g
Chromogenic solution	1 ml
Selective solution	20 ml
Agar B	13 g
Distilled water	qsp 1,000 ml

---

Final pH at 25°C = 7.2 ± 0.2

## Section 4 Shelf Life and Storage

- Dehydrated agar: 2–8°C in carefully sealed package, in a dry and dark place
- Pre-poured agar: 2–8°C in carefully sealed package, in a dark place
- Plate prepared from the dehydrated agar: 1 week at 2–8°C in carefully sealed package, in a dry and dark place

## Section 5 Materials Required but Not Supplied

### Equipment

- All usual laboratory equipment
- Hot plate
- Scale, sensitivity of 0.1 g
- Stirrer/homogenizer
- Thermostatically controlled incubator or incubation room, precise to  $\pm 1^\circ\text{C}$
- Water bath

### Supplies

- Confirmation:
  - iQ-Check *Listeria monocytogenes* II Real-Time PCR Detection Kit (catalog #3578124), iQ-Check *Listeria* spp. Real-Time PCR Detection Kit (catalog #3578113)
  - Agar *Listeria* according to Ottaviani and Agosti (for example, catalog #3563695, prepared plates 90 mm x 20; 3563965, prepared plates 90 mm x 120; 3555200, 250 ml x 6 bottles; 3564043, dehydrated, 500 g; 3564041, AL Supplement 1, 10 vials; 3564042, AL Supplement 2, 10 vials)
  - PALCAM agar (for example catalog #3564754, dehydrated 500 g; 3564752, supplement, 10 vials)
  - Rhamnose test (catalog #3553669, 1 ml x 28 vials)
- Diluent for enumeration, tryptone salt (catalog #3555754, 9 ml x 25 tubes; 3555756, 900 ml x 6 bottles; 3555796, 3 L x 4 bags; 3564544, 500 g)
- Enrichment medium, Fraser ½ broth (catalog #3555797, 225 ml bottles x 6; 3555794, 3 L x 4 bags; 3564604, dehydrated 500 g; 3564616, supplement, 10 vials)
- Inoculating loops and spreaders

## Section 6 Precautions, Limitations of Use, and Quality Control

- Resuscitation medium, Buffered Peptone Water (BPW), (catalog #3554179, 225 ml x 6 bottles; 3564684, dehydrated 500 g; 3555790, 2 x 5 L bags; 3555795, 4 x 3 L bags)
- Sterile petri dishes
- Sterile pipets
- Sterile weigh bags

## Section 6 Precautions, Limitations of Use, and Quality Control

### Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria such as *L. monocytogenes*
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Avoid prolonged overheating of the product during melting
- Before using RAPID'*L.mono* Agar, let the petri dishes dry at 25–50°C until droplets disappear from the surface of the medium, in compliance with the EN ISO 7218 standard. Avoid prolonged drying, however, as this could alter the performance of the medium
- When using the enumeration method, inoculate the medium using a spreader. After spreading, the dishes can be left as they are on the work surface for 15–30 min before incubation in order to permit the inoculum to be correctly absorbed by the agar
- When using bottled RAPID'*L.mono* Agar, a white deposit at the bottom of the bottle is normal. To ensure resuspension and satisfactory homogeneity, it is important to stir the bottle by hand when it is taken out of the boiling and melting water baths and to pour it immediately after adding the supplements
- Before inoculation, RAPID'*L.mono* Agar is red to red-orange and opaque. Some color variation can be observed without any impact on the performance

### Limitations of Use

- A strain of *L. ivanovii* with slow xylose activity can be found in sheep milk. This typical strain is difficult to differentiate from a *L. monocytogenes* strain, even after 48 hr on RAPID'*L.mono* Agar. Therefore, Agar *Listeria* according to Ottaviani and Agosti confirmation isn't advised when sheep milk products are tested
- Some rare *L. monocytogenes* strains are characterized by slow PIPLC activity, which could result in a slower development of the typical blue color. Extending the incubation time by 48 hr may be

## Section 7 Protocol

required. During the NF Validation study, 200 strains of *L. monocytogenes* were tested and none demonstrated slow PIPLC activity

- Sample volumes over 25 g have not been tested as part of NF VALIDATION

## Quality Control

- Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file
- For SDS product safety information and certificate of analysis, visit [www.bio-rad.com](http://www.bio-rad.com)

## Section 7 Protocol

### Preparation of Dehydrated Medium

1. Always shake bottle before use.
2. Dissolve 68.1 g of powder in 950 ml of distilled water and mix until a homogenous suspension is obtained.
3. Warm slowly while stirring and bring to a boil (maximum boil 2 min). Avoid overheating.
4. Cool the medium to 45–50°C in a water bath.
5. Aseptically add and homogenize 2 vials of Supplement 1, each reconstituted by adding 14 ml of sterile water.
6. Aseptically add and homogenize 2 vials of Supplement 2.
7. If necessary, mix on a magnetic stirrer until dispensed.
8. Dispense in petri dishes.
9. Leave to set on a cool, level surface.
10. Do not stack the dishes.
11. One 500 g bottle of powder makes 7.4 L of medium.

### Preparation of Bottled Medium

1. In a boiling water bath, melt the contents of the RAPID'*L.mono* Agar bottle.

## Section 7 Protocol

2. Remove from heat and stir the bottle contents by hand to resuspend the white deposit.
3. Cool the bottle to 47–50°C in a water bath. Stir bottle contents by hand before adding supplements.
4. Aseptically add 1 vial of Supplement 1 reconstituted with 14 ml of sterile water and 1 vial of Supplement 2.
5. Mix thoroughly. Avoid frothing.
6. Immediately pour the medium into petri dishes.
7. Leave to set on a cool, level surface.
8. Do not stack the dishes.
9. One kit allows the preparation of approximately 13 plates.

### **Detection of *L. monocytogenes* and *Listeria* genus**

1. Dilute  $n$  g or  $n$  ml of sample in  $9 \times n$  ml Fraser  $\frac{1}{2}$  broth. Homogenize with stirrer/homogenizer.
2. Incubate at  $30 \pm 1^\circ\text{C}$  for  $24 \pm 2$  hr.
3. After the enrichment step, the broth can be stored at  $2\text{--}8^\circ\text{C}$  for 72 hr.
4. Using a sterile pipet remove 0.1 ml of sample and place drops on the outside edge of half of the agar surface.
5. Using a sterile Pasteur pipet or inoculating loop, spread sample over half the agar surface in a to-and-fro motion.
6. On the other half of the agar surface, streak for isolation by spreading the deposit in relatively close streaks over the entire dish from the edge of the previous spread.
7. Incubate the upturned dishes at  $37 \pm 1^\circ\text{C}$  for  $24 \pm 2$  hr.
8. *L. monocytogenes* form pale blue, grey-blue to dark blue round and smooth colonies, without yellow halo, with an average diameter of 1–2 mm. Depending on the food matrix, the depth of blue color of the colony may vary.
9. Other *Listeria* spp. colonies are typically white or pale yellow, with or without a yellow halo, forming a round shape with a smooth, convex appearance, average diameter 1–2 mm.  
*L. ivanovii* form blue-green colonies with a yellow halo.
10. In the case of "non-detected" target of *Listeria* genus (other than *L. monocytogenes*), 48 hr of incubation are required to provide a negative result.
11. After the incubation step, plates can be stored at  $2\text{--}8^\circ\text{C}$  for 72 hr.

## Section 8 Confirmation of Positive Results

### Enumeration of *L. monocytogenes*

1. Dilute  $n$  g or ml of sample in  $9 \times n$  ml of BPW.
2. Incubate at  $20 \pm 1^\circ\text{C}$  for  $1 \text{ hr} \pm 5 \text{ min}$ . This step is optional; the resuscitation step was not applied during the 2018 NF VALIDATION renewal study.
3. Spread 0.1 ml on one plate.
4. If necessary, prepare a 1:10 dilution (or more) in tryptone salt diluent or BPW as per the ISO 6887 standard and spread 0.1 ml of each dilution on a plate.
5. If estimating small numbers, spread 1 ml of sample over three plates of 90 mm diameter ( $\sim 0.33$  ml/plate) or over one plate of 140 mm diameter (NF VALIDATION certified protocol). It is also possible to spread 1 ml over two plates of 90 mm diameter (not NF VALIDATION certified).
6. Incubate the upturned dishes at  $37 \pm 1^\circ\text{C}$  for  $24 \pm 2 \text{ hr}$ .
7. Read plate and enumerate typical colonies. Count only dishes containing a maximum of 150 characteristic colonies and a maximum of 300 colonies in total.
8. Refer to EN ISO 7218 standard for inoculation, colony counting, calculation and expression of results.

## Section 8 Confirmation of Positive Results

1. In the context of AOAC validation, confirm suspect isolated colonies according to classic confirmation test procedure described in the standard reference methods.
2. In the context of NF VALIDATION, confirmation of *L. monocytogenes* performed on only one colony is sufficient and must be done in one of the following ways:
  - a. Using the conventional tests described in the CEN or ISO standard methods (including the purification step).
  - b. Using nucleic probes as described in the ISO 7218 standard (for example, iQ-Check *L. monocytogenes* II Real-Time PCR Detection Kit, catalog #3578124) using isolated colonies (with or without purification step).
  - c. Confirmed using a rhamnose test (catalog # 3553669).
  - d. Using any other NF VALIDATION certified method based on a different principle from that of RAPID'*L.mono*. The validated protocol of the second method must be respected in its entirety. All steps preceding the detection step used as a starting point for confirmation must be common to both methods.
3. In the context of NF VALIDATION, confirmation of *Listeria* spp. other than *L. monocytogenes* performed on only one colony is sufficient and must be done in one of the following ways:

## Section 9 Confirmation of Other Methods

- a. Using the conventional tests described in the CEN or ISO standard methods (including the purification step).
  - b. Using nucleic probes as described in the ISO 7218 standard (for example, iQ-Check *Listeria* spp. Real-Time PCR Detection Kit, catalog #3578113) using isolated colonies (with or without purification step).
  - c. Repicking and spotting at least one isolated colony on an AL or a PALCAM agar plate. Up to 15 colonies can be confirmed on a single AL or PALCAM agar plate.
  - d. Using any other NF VALIDATION certified method based on a different principle from that of RAPID'*L.mono* Agar. The validated protocol of the second method must be respected in its entirety. All steps preceding the detection step used as a starting point for confirmation must be common to both methods.
4. Suspected colonies can be confirmed by using the MALDI Biotyper test (Bruker's Biotyper System includes a MALDI time-of-flight mass spectrometer and MBT Compass Software, version 4) directly from an isolated colony or after a purification step. This option is certified by MicroVal according to ISO 16140-6 (2017LR75).
  5. In the event of discordant results (presumptive positive with RAPID'*L.mono*, negative with confirmation method), the laboratory must follow the necessary steps to ensure validity of the result obtained.
  6. In the case of a positive result confirmed for *L. monocytogenes* with the RAPID'*L.mono* enumeration method, it is not necessary to duplicate the *L. monocytogenes* confirmation step for the same sample tested with the RAPID'*L.mono* detection method.
  7. When using the enumeration protocol, confirming fewer than five colonies involves a risk of overestimating because of the presence of typical colonies that would not be *L. monocytogenes*.

## Section 9 Confirmation of Other Methods

### Confirmation of Positive Results with iQ-Check *Listeria* spp. and iQ-Check *Listeria monocytogenes* II Kits (NF VALIDATION Certified Protocol)

1. Using a sterile loop, streak 100 µl of enrichment broth (LSB or Fraser ½ broth) on RAPID'*L.mono* Agar.
2. Incubate plates at 37 ± 1°C for 24 ± 2 hr.
3. *L. monocytogenes* form pale blue, grey-blue to dark blue colonies and other *Listeria* spp. form white or pale yellow colonies with or without a yellow halo on RAPID'*L.mono* Agar.

## Section 10 Test Performance and Validations

Certification	Scope	Validation Protocol	Reference Protocol	Certificate Reference
AOAC-RI	Brie cheese, surimi, mixed salad, deli turkey	Performance Tested Methods	FDA BAM Chap. 10 USDA MLG 8.11	 License# 030406
NORDVAL	All human food products and production environmental samples	EN ISO 16140-2	EN ISO 11290-1 (2017) EN ISO 11290-2 (2017)	 NordVal #022
NF VALIDATION	All human food products and production environmental samples	EN ISO 16140-2	EN ISO 11290-1 (2017) EN ISO 11290-2 (2017)	 BRD 07/04–09/98 BRD 07/05–09/01 ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS Certified by AFNOR Certification <a href="http://nf-validation.afnor.org/en">http://nf-validation.afnor.org/en</a> Refer to the certificates available on the website mentioned above for the end of validity.

## Section 11 References

Amilli A, Goldfine H, and Portnoy DA. (1991). *Listeria monocytogenes* mutants lacking phosphatidylinositol-specific phospholipase C are avirulent. *J Exp Med* 173, 751–754.

Becker B, Schuler S, Lohneis M, Sabrowski A, Curtis D.W, Holzapfel W. (2006). Comparison of two chromogenic media for the detection of *Listeria monocytogenes* with the plating media recommended by EN/DIN 11290-1. *Int J Food Microbiol* 109, 127–131.

Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. (1990). Enhanced haemolysis agar (EHA)—an improved selective and differential medium for isolation of *Listeria monocytogenes* *Food Microbiol* 8, 37–49.

## Section 11 References

Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. 1990: Performance of enhanced haemolysis agar compared to Oxford medium for the isolation of *Listeria monocytogenes* from food enrichment broths. *Food Microbiol* 8, 51-62.

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Heisick JE, Rosas-Marty LI and Tatini SR. (1995). Enumeration of Viable *Listeria* Species and *Listeria monocytogenes* in Foods. *J Food Prot* 58, 733–736.

ISO 11290-1:2017 - Microbiology of the food chain -- Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. -- Part 1: Detection method.

ISO 11290-2:2017 - Microbiology of the food chain -- Horizontal method for the detection and enumeration of *Listeria monocytogenes* and of *Listeria* spp. -- Part 2: Enumeration method.

ISO 6887-1:2017/ADM1:2024. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: General rules for the preparation of the initial suspension and decimal dilutions

ISO 18593: 2018. Microbiology of the food chain - Horizontal methods for surface sampling

Karpiskova R., Pejchalova M., Mokrosova J, Vytrasova J, Smuharova P, Ruprich J. (2000). Application of a chromogenic medium and the PCR method for the rapid confirmation of *L. monocytogenes* in foodstuffs. *J Microbiol Methods* 41, 267-271.

Leclercq A (2004). Atypical colonial morphology and low recoveries of *Listeria monocytogenes* strains on Oxford, PALCAM, Rapid'*L.mono* and ALOA solid media. *J Microbiol Methods* 57, 251–258.

Manafi M (2000). New developments in chromogenic and fluorogenic culture media. *Int J Food Microbiol* 60, 205-218.

Mengaud J, Braun-Breton C, and Cossart P. (1991). Identification of phosphatidylinositol-specific phospholipase C activity in *Listeria monocytogenes*: a novel type of virulence factor? *Mol Microbiol*, 5, 367–372.

Notermans SH, Dufrenne J. et al. (1991). Phosphatidylinositol-specific phospholipase C activity as a marker to distinguish between pathogenic and nonpathogenic *Listeria* species. *Appl Envir Microbiol* 57, 2666–2670.

Polivka C (2001). Identification of *Listeria* with a new chromogenic medium RAPID'*L.mono*. *Arch Lebensmittelhygiene* 52, 22–23.

United States Department of Agriculture. *Microbiology Laboratory Guidebook*. Chapter 8.11 Isolation and Identification of *Listeria monocytogenes* from Red Meat, Poultry, Ready-to-Eat Siluriformes (Fish) and Egg Products, and Environmental Sponges. Online at <https://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/mlg-8.pdf?MOD=AJPERES>

United States Food and Drug Administration. *Bacteriological Analytical Manual*. Chapter 10 *Listeria monocytogenes*. Online at <https://www.fda.gov/food/laboratory-methods-food/bacteriological-analytical-manual-bam>

## Section 12 Revision History

Release date	Document number	Change
March 2025	10000251398 Ver A	- Document number change (previous version 10000127436 Ver C) - References update

Visit [www.bio-rad.com/rapidmedia](http://www.bio-rad.com/rapidmedia) for more information on our complete range of RAPID Chromogenic Media.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. IQ-CHECK is a trademark of Bio-Rad Europe GmbH in certain jurisdictions. All trademarks used herein are the property of their respective owner.

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Becker B, Schuler S, Lohneis M, Sabrowski A, Curtis D.W, Holzapfel W. (2006). Comparison of two chromogenic media for the detection of *Listeria monocytogenes* with the plating media recommended by EN/DIN 11290-1. *Int J Food Microbiol* 109, 127-131.

Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. (1990). Enhanced haemolysis agar (EHA)—an improved selective and differential medium for isolation of *Listeria monocytogenes* *Food Microbiol* 8, 37-49.

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Heisick JE, Rosas-Marty LI and Tatini SR. (1995). Enumeration of Viable *Listeria* Species and *Listeria monocytogenes* in Foods. *J Food Prot* 58, 733-736.

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Leclercq A (2004). Atypical colonial morphology and low recoveries of *Listeria monocytogenes* strains on Oxford, PALCAM, Rapid'Lmono and ALOA solid media. *J Microbiol Methods* 57, 251-258.

Manafi M (2000). New developments in chromogenic and fluorogenic culture media. *Int J Food Microbiol* 60,205-218.

Mengaud J, Braun-Breton C, and Cossart P. (1991). Identification of phosphatidylinositol-specific phospholipase C activity in *Listeria monocytogenes*: a novel type of virulence factor? *Mol Microbiol*, 5, 367-372.

Notermans SH, Dufrenne J. et al. (1991). Phosphatidylinositol-specific phospholipase C activity as a marker to distinguish between pathogenic and nonpathogenic *Listeria* species. *Appl Environ Microbiol* 57, 2666-2670.

Polivka C (2001). Identification of *Listeria* with a new chromogenic medium RAPID'/.*mono*. *Arch Lebensmittelhygiene* 52, 22-23.

United States Department of Agriculture. *Microbiology Laboratory Guidebook*. Chapter 8.11 Isolation and Identification of *Listeria monocytogenes* horn Red Meat, Poultry, Ready-to-Eat Siluriformes (Fish) and Egg Products, and Environmental Sponges. Online at <https://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/mlq-8.pdf?MOD=AJPERES>

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Becker B, Schuler S, Lohneis M, Sabrowski A, Curtis D.W, Holzapfel W. (2006). Comparison of two chromogenic media for the detection of *Listeria monocytogenes* with the plating media recommended by EN/DIN 11290-1. Int J Food Microbiol 109, 127-131.

Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. (1990). Enhanced haemolysis agar (EHA)—an improved selective and differential medium for isolation of *Listeria monocytogenes* Food Microbiol 8, 37-49.

Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. 1990: Performance of enhanced haemolysis agar compared to Oxford medium for the isolation of *Listeria monocytogenes* from food enrichment broths. Food Microbiol 8, 51-62.

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Heisick JE, Rosas-Marty LI and Tatini SR. (1995). Enumeration of Viable *Listeria* Species and *Listeria monocytogenes* in Foods. J Food Prot 58, 733-736.

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Manafi M (2000). New developments in chromogenic and fluorogenic culture media. Int J Food Microbiol 60,205-218.

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Notermans SH, Dufrenne J. et al. (1991). Phosphatidylinositol-specific phospholipase C activity as a marker to distinguish between pathogenic and nonpathogenic *Listeria* species. *Appl Envir Microbiol* 57, 2666-2670.

Polivka C (2001). Identification of *Listeria* with a new chromogenic medium RAPID'*L. mono*. *Arch Lebensmittelhygiene* 52, 22-23.

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ISO 6887-1:2017/ADM1:2024. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: General rules for the preparation of the initial suspension and decimal dilutions.

ISO 18593: 2018. Microbiology of the food chain - Horizontal methods for surface sampling.

Karpiskova R., Pejchalova M., Mokrosova J, Vytrasova J, Smuharova P, Ruprich J. (2000). Application of a chromogenic medium and the PCR method for the rapid confirmation of *L. monocytogenes* in foodstuffs. *J Microbiol Methods* 41, 267-271.

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United States Department of Agriculture. *Microbiology Laboratory Guidebook*. Chapter 8.11 Isolation and Identification of *Listeria monocytogenes* horn Red Meat, Poultry, Ready-to-Eat Siluriformes (Fish) and Egg Products, and Environmental Sponges. Online at <https://www.fsis.usda.gov/wps/wcm/connect/1710bee8-76b9-4e6c-92fc-fdc290dbfa92/mlg-8.pdf?MOD=AJPERES>

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Cox LJ, Siebenga A, Pedrazzini C, and Moreton J. 1990: Performance of enhanced haemolysis agar compared to Oxford medium for the isolation of *Listeria monocytogenes* from food enrichment broths. *Food Microbiol* 8, 51–62.

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ISO 6887-1:2017. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: General rules for the preparation of the initial suspension and decimal dilutions

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Manafi M (2000). New developments in chromogenic and fluorogenic culture media. *Int J Food Microbiol* 60, 205–218.

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