
RAPID'*Enterobacteriaceae* Agar

User Guide

Chromogenic media for the detection and enumeration of *Enterobacteriaceae*, without confirmation, in 24 hr in human food, feed, and environmental samples

Catalog #3554012, Ready-to-use, 200 ml x 6 bottles
Catalog #3564004, Dehydrated, 500 g

BIO-RAD

Table of Contents

Section 1	Introduction	1
Section 2	RAPID' <i>Enterobacteriaceae</i> Principle.....	1
Section 3	Theoretical Formula.....	1
Section 4	Shelf Life and Storage.....	1
Section 5	Materials Required but Not Supplied	2
	Equipment	2
	Supplies	2
Section 6	Precautions, Limitations of Use, and Quality Control	2
Section 7	Protocol	3
	Preparation of Dehydrated Medium.....	3
	Sample Preparation	4
	Inoculation and Plate Reading	4
Section 8	Confirmation of Positive Results	4
Section 9	Confirmation of Other Methods.....	4
Section 10	Test Performance and Validations.....	5
Section 11	References	5
Section 12	Revision History	5

Section 1 Introduction

Indicator microorganisms are used as a sign of quality or hygienic status in food, feed, water, and environmental samples. These microorganisms can signify the presence of pathogens and highlight issues with sanitation or production processes. Different indicators may be used for different applications, according to context and goals. *Enterobacteriaceae* is a large family of rod-shaped gram-negative bacteria. Some are commensal; some, such as *Salmonella*, *Cronobacter*, and some *Escherichia coli* strains are pathogenic. Quantification of *Enterobacteriaceae* is often used for hygiene and quality monitoring in food matrices and environmental samples.

Section 2 RAPID'*Enterobacteriaceae* Principle

The principle of the RAPID'*Enterobacteriaceae* medium relies on the ability of *Enterobacteriaceae* to ferment glucose. Due to the simultaneous presence of crystal violet and bile salts, the medium inhibits gram-positive bacteria and some gram-negative bacteria. The combination of color indicators allows a high level of contrast of *Enterobacteriaceae* colonies which appear red on a clear gray medium.

Section 3 Theoretical Formula

Nutritive mix	17.3 g
Glucose	9.0 g
Selective agents	0.7 g
Color indicators	0.07 g
Agar	11 g
Distilled water	qsp 1,000 ml

Final pH at 25°C = 7.4 ± 0.2

Section 4 Shelf Life and Storage

- Dehydrated: 15–25°C in carefully sealed package, in a dry and dark place
- Bottles: 2–8°C in a dark place
- Plate prepared from the dehydrated agar: 2 weeks at 2–8°C in carefully sealed package, in a dry and dark place

Section 5

Materials Required but Not Supplied

Equipment

- All usual laboratory equipment
- Hot plate
- Scale, sensitivity of 0.1 g
- Stirrer/homogenizer
- Thermostatically-controlled incubator or incubation room, precise to $\pm 1^{\circ}\text{C}$
- Water bath

Supplies

- Diluent for enumeration: Buffered Peptone Water (BPW) (catalog #3554179, 225 ml x 6 bottles; 3564684, dehydrated, 500 g; 3555790, 2 x 5 L bags; 3555795, 4 x 3 L bags); Tryptone Salt (catalog #3555754, 9 ml x 25 tubes; 3555756, 90 ml x 6 bottles; 3564544, dehydrated, 500 g; 3555796, 3 L x 4 bags)
- Inoculating loops
- Sterile petri dishes (90 mm)
- Sterile pipets
- Sterile weigh bags

Section 6

Precautions, Limitations of Use, and Quality Control

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria such as *Enterobacteriaceae*
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations
- Do not autoclave the medium

Section 7 Protocol

- It is preferable to use a double-layered medium in matrices containing abundant mesophilic flora. The aim of the second layer is to limit invasion of the surface, which can interfere with reading

Limitations of Use

- In the context of the NF VALIDATION mark, a spiral plate system was not been tested
- In the context of the NF VALIDATION mark, Scan 1200 instrument with software version "V 6.3.7 kernel Biorad 1.2" was tested. Refer to the instrument user manual for any questions regarding the use of the colony counter. A median bias of 0.23 log CFU/g was observed during the validation study using the surface inoculation method with the automatic colony counter

Quality Control

- Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file.
- For SDS product safety information and certificate of analysis, visit www.bio-rad.com.

Section 7 Protocol

Preparation of Dehydrated Medium

1. Always shake bottle before use.
2. Dissolve 38 g of powder in 1 L of distilled water and mix until a homogenous suspension is obtained.
3. Heat gently, agitating frequently, then bring to a boil.
4. Do not autoclave.
5. Cool the medium to 44–47°C before use.
6. This medium can be used directly after preparation or dispensed in bottles that are sealed and stored in a dry cool place. Sterility of the medium should be checked before use.
7. One 500 g bottle of powder makes 13.5 L of medium.

Sample Preparation

Dilute sample according to the standard method applicable to the product concerned. In the context of the NF Validation certified method, prepare the initial suspension and subsequent dilutions in accordance with the ISO 6887 standard series.

Inoculation and Plate Reading

1. **Pour plate method:** Pour 1 ml of sample and/or its decimal dilutions in single layer in an empty petri dish. Pour the melted medium, cooled to 44–47 °C, onto the sample and homogenize by swirling. A second layer can be poured (approximately 2 mm thick) in case of high level of contamination.
2. **Surface inoculation method:** Spread 0.1 ml of sample and/or its decimal dilutions on a petri dish. If it is necessary to estimate small numbers, spread 1 ml of sample over three dishes of Ø 90 mm (~ 0.33 ml/dish) or over one dish of 140 mm.
3. Incubate plates at 37 ± 1°C for 24 ± 2 hr. Alternatively, a temperature of 30 or 35°C can be used when enumerating specific *Enterobacteriaceae* (e.g., psychrotrophic *Enterobacteriaceae*).
4. *Enterobacteriaceae* form red colonies with a diameter equal to or exceeding 0.5 mm, with or without a zone of precipitation.
5. **Pour plate method:** The reading step can be also performed with the Scan 1200 Automatic Colony Counter with the pour plate inoculation protocol. Select "Rapid Entero" and follow the instructions in order to enumerate the colonies.
6. Refer to ISO 7218 standard for inoculation, colony counting, calculation and expression of results. Enumeration using a single plate and a single dilution is an option certified under NF Validation by AFNOR Certification. In this context, the requirements of ISO 7218 concerning the use of two successive dilutions or two plates from the same dilution cannot be applied (see the sections "Inoculation" and "Calculations and Expression of Results").
7. **Pour plate method:** In the context of the NF Validation certified method, after the incubation step, plates can be stored at 2–8°C for 72 hr prior to colony counting step.
8. **Surface inoculation method:** After the incubation step, plates can be stored at 2–8°C for 48 hr according to ISO 7218

Section 8 Confirmation of Positive Results

Not applicable.

Section 9 Confirmation of Other Methods

Not applicable.

Section 10 Test Performance and Validations

Certification	Scope	Validation Protocol	Reference Protocol	Certificate Reference
NF VALIDATION	All human food, feed, and environmental samples	EN ISO 16140-2	ISO 21528-2	 <p>BRD 07/24-11/13 ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS Certified by AFNOR Certification http://nf-validation.afnor.org/en Refer to the certificate available on the website mentioned above for the end of validity.</p>

Section 11 References

ISO 6887-1:2017 /AMD1: 2024. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 1: General rules for the preparation of the initial suspension and decimal dilutions

ISO 7218:2024. Microbiology of the food chain – General requirements and guidance for microbiological examinations

ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 2: Colony-count technique

Section 12 Revision History

Release date	Document number	Change
June 2025	10000251395 Ver A	<ul style="list-style-type: none"> - Document number change (previous version 10000127900 Ver B) - Extension of NF Validation for simple plate and single dilution protocol - Update of bibliographic references

Visit www.bio-rad.com/rapidmedia for more information on our complete range of RAPID chromogenic media.

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

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BIO-RAD es una marca registrada de Bio-Rad Laboratories, Inc.