RAPID'*B.cereus* Agar

User Guide


Chromogenic media for the detection and the enumeration of *Bacillus cereus* group (vegetative cells and spores) in human food products, animal feed and environmental samples at 30°C in 24 hr without confirmation

Catalog #12007304, Prepared plates, 90 mm x 20 dishes

Catalog #12007305, Dehydrated, 500 g

Catalog #12007306, Supplement 1, liquid, 10 vials

Catalog #12007307, Supplement 2, freeze-dried, 10 vials

BIO-RAD

Table of Contents

Section 1	Introduction	1
Section 2	RAPID' <i>B.cereus</i> Principle	1
Section 3	Theoretical Formula.....	1
Section 4	Shelf Life and Storage.....	1
Section 5	Materials Required but Not Supplied	2
	Equipment	2
	Supplies	2
Section 6	Precautions, Limitations of Use, and Quality Control	2
Section 7	Protocol	3
	Preparation of Dehydrated Medium.....	3
	Sample Preparation	4
	Inoculation and Plate Reading	4
Section 8	Confirmation of Positive Results	4
Section 9	Confirmation of Other Methods.....	4
Section 10	Test Performance and Validations.....	5
Section 11	References	5
Section 12	Revision History	5

Section 1 Introduction

The *Bacillus cereus* group of organisms is widely distributed in the environment and eliminating it can present a challenge in food production, particularly with ready-to-eat and ready-to-reheat products. Different bacteriological characteristics, including biochemical characteristics, are shared by several species within the *Bacillus* genus, making differentiation of these organisms difficult. These spore-forming organisms have traditionally been controlled by the application of heat. However, the increasing demand for food with a fresh-like quality (texture, color, and aroma) and the desire to produce foods with functional ingredients have increased interest in the development of low-heat operations to either inactivate spores or to prevent their germination and growth.

Section 2 RAPID'*B.cereus* Principle

RAPID'*B.cereus* Agar is based on a chromogenic reaction and phospholipase activity. Typical *B. cereus* colonies develop a characteristic red color generally surrounded with an opaque halo. Some rare *B. cereus* do not express the phospholipase activity and aren't able to produce the opaque halo. The medium's selective mixture prevents the growth of the interfering flora and enables the analysis of a broad range of foods. The sensitivity of the medium has been specially optimized to allow growth of all pathogenic *Bacillus cereus* strains, even those that are difficult to culture, like *B. cytotoxicus*. The interpretation of RAPID'*B.cereus* Agar plates is further facilitated by the high level of contrast between the color of the colonies and the agar.

Section 3 Theoretical Formula

Peptones	13 g
Phospholipids	5 g
Selective agents	1.1 g
Chromogenic mix	0.01 g
Salt	10.5 g
Agar	15 g
Distilled water	qsp 1,000 ml

Final pH at 25°C = 7.0 ± 0.2

Section 4 Shelf Life and Storage

- Dehydrated agar: 15–25°C in a carefully sealed package in a dry and dark place
- Supplements: 2–8°C in a dark place

Section 5 Materials Required but Not Supplied

- Preprepared agar plates: 2–8°C in a dark place
- Plate prepared from the dehydrated agar: 1 month at 2–8°C in a carefully sealed package in a dry and dark place

Section 5 Materials Required but Not Supplied

Equipment

- All usual laboratory equipment
- Hot plate
- Scale, sensitivity of 0.1 g
- Stirrer/homogenizer
- Thermostatically controlled incubator or incubation room, precise to $\pm 1^\circ\text{C}$
- Water bath

Supplies

- Diluent for enumeration: Buffered Peptone Water (BPW) (e.g. catalog #3554179, 225 ml x 6 bottles; 3564684, dehydrated, 500 g; 3555790, 5 L x 2 bags; 3555795, 3 L x 4 bags); tryptone salt (e.g. catalog #3555754, 9 ml x 25 tubes; 3555756, 90 ml x 6 bottles; catalog #3564544, dehydrated, 500 g; catalog #3555796, 3 L x 4 bags)
- Inoculating loops
- Sterile petri dishes
- Sterile pipets
- Sterile weigh bags

Section 6 Precautions, Limitations of Use, and Quality Control

Precautions

- Respect Good Laboratory Practice (EN ISO 7218). Appropriate protection, such as gloves and lab coats, should be worn when working with potentially infectious live bacteria such as *Bacillus*
- Media that have come in contact with food samples should be considered contaminated and should be disposed of in accordance with local rules and regulations

Section 7 Protocol

- The time lapse between the end of preparation of the stock solution (or 10⁻¹ dilution in the case of a solid product) and when the dilutions come into contact with the culture medium must not exceed 45 min
- Plates should be pre-dried before use to ensure that they are free of excess moisture (EN ISO 7218). For example, plates can be dried by placing them with the agar surface facing upwards in a laminar-flow safety cabinet (at room temperature) for 20 to 30 minutes.

Limitations of Use

- A small percentage of *Bacillus cereus* strains do not express phospholipase activity (red colonies without opaque halo), especially when the pour plate protocol is used. These red colonies should be considered typical

Quality Control

- Every product manufactured and marketed by Bio-Rad is subject to a quality assurance procedure at all stages, from reception of raw materials through to marketing of the finished products. Each batch of finished product undergoes quality control according to EN ISO 11133 and is marketed only if it satisfies the acceptability criteria. Documentation relative to the production and quality control of each batch is kept on file
- For SDS product safety information and certificate of analysis, visit www.bio-rad.com

Section 7 Protocol

Preparation of Dehydrated Medium

1. Always shake bottle before use.
2. Dissolve 36.5 g of powder in 950 ml of distilled water, wait 5 min, then mix until a homogenous suspension is obtained.
3. Heat gently, agitating frequently, then bring to a boil.
4. Dispense 475 ml into two bottles.
5. Sterilize by autoclaving at 121°C for 15 min.
6. Cool the medium to 44–47°C.
7. Aseptically add 1 vial of Supplement 1.
8. Aseptically add 1 vial of Supplement 2 that has been reconstituted in 6 ml sterile water.
9. Mix well and pour into petri dishes or keep in bottles to use a pour plate inoculation method.
10. Complete medium prepared with supplements can be maintained at 44–47°C up to 8 hr before use.
11. One 500 g bottle of powder makes 13.8 L of medium.

Section 8 Confirmation of Positive Results

Sample Preparation

Prepare sample according to the standard applicable to the product concerned.

In the context of an NF Validation certified method, prepare the initial suspension and subsequent dilutions in accordance with ISO 6887-1 and the relevant parts of ISO 6887 that apply to the specific sample type.

Inoculation and Plate Reading

1. **Pour plate method:** Transfer 1 ml of sample and/or its decimal dilutions into an empty petri dish. Pour the melted medium, cooled to 44–47 °C, onto the sample and homogenize by swirling. Let solidify on a cool surface.
2. **Surface inoculation method:** Spread 0.1 ml of sample and/or its decimal dilutions on a petri dish. If it is necessary to estimate small numbers, spread 1 ml of sample over three 90 mm dishes (~0.33 ml/dish) or over one 140 mm dish.
3. For both methods, incubate plates facing down at $30 \pm 1^\circ\text{C}$ for 24 ± 3 hr.
4. *B. cereus* produces red colonies surrounded with an opaque halo or without halo in some rare cases. The aspect of strains belonging to the *B. cereus* group may be variable (rhizoidal aspect, rough colonies, variable diameter, etc.).
5. Refer to EN ISO 7218 standard for inoculation, colony counting, calculation and expression of results. Enumeration using a single plate and a single dilution is an option certified under NF Validation by AFNOR Certification. In this context, the requirements of ISO 7218 concerning the use of two successive dilutions or two plates from the same dilution cannot be applied (see the sections "Inoculation" and "Calculations and Expression of Results").
6. After the incubation step, plates can be stored at 2–8°C for 72 hr for surface inoculation.

Section 8 Confirmation of Positive Results

Not applicable.

Section 9 Confirmation of Other Methods

Not applicable.

Section 10 Test Performance and Validations

Certification	Scope	Validation Protocol	Reference Protocol	Certificate Reference
NF VALIDATION	All food products, animal feed, environmental samples	EN ISO 16140-2	ISO 7932	 <p>BRD 07/26-03/19 ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS Certified by AFNOR Certification http://nf-validation.afnor.org/en Refer to the certificate available on the website mentioned above for the end of validity.</p>

Section 11 References

ISO 7932:2004/AMD 1:2020 – Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of presumptive *Bacillus cereus* – Colony-count technique at 30 °C.

ISO 7218:2024. Microbiology of the food chain – General requirements and guidance for microbiological examinations

ISO 6887-1:2017 /AMD1: 2024. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination – Part 1: General rules for the preparation of the initial suspension and decimal dilutions

Section 12 Revision History

Release date	Document number	Change
June 2025	10000251397 Ver A	<ul style="list-style-type: none"> - Document number change (previous version 10000127897 Ver C) - Extension of NF Validation for simple plate and single dilution protocol - Update of bibliographic references

Visit www.bio-rad.com/rapidmedia for more information on our complete range of RAPID chromogenic media

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.

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ISO 7932:2004/AMD 1:2020 — Microbiology of food and animal feeding stuffs — Horizontal method for the enumeration of presumptive *Bacillus cereus* — Colony-count technique at 30 °C.

ISO 7218:2024. Microbiology of the food chain - General requirements and guidance for microbiological examinations

ISO 6887-1:2017 /AMD1: 2024. Microbiology of the food chain - Preparation of test samples, initial suspension and decimal dilutions for microbiological examination - Part 1: General rules for the preparation of the initial suspension and decimal dilutions

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