

Profinia User Interface (UI) Firmware Upgrade Instructions

To download Profinia UI firmware, version 1.40 to your Profinia instrument:

1. Turn your computer and Profinia instrument on.
2. Verify that the USB cable is connected from the Profinia instrument to your computer.



3. Close all open Windows applications on the computer, including the Profinia software, before starting the firmware download.
4. Locate the file “Profinia UI Firmare Ver 140.exe” on your computer.



5. Double-click the file to open the Profinia User Interface Firmware Ver. 1.40 dialog box.



6. Click **Download Firmware**. The firmware download and verification will take approximately 18 minutes. The Progress field shows percentage complete.
7. When the message “Verify Finished. Cycle Power” is displayed, the firmware download and the verification is finished. Close the program and turn off the Profinia instrument. Your Profinia instrument is now ready to use.

Important note regarding Saved Program Methods:

The **Saved Methods** created using older versions of Profinia firmware will not work properly using the new firmware version 1.40. We recommend that you record the parameters saved in the old methods and reprogram them as new methods on the new firmware.

To print out the steps of a Saved Method using Profinia software:

The run files that have been saved in Profinia software will not be affected by the firmware change. To print out the detailed steps of a Saved Method, use Profinia software to open a purification run that used the method you want to reprogram. To print out the steps:

- Select “Reports” from the upper toolbar
- Select “Custom Reports”
- Select the “Method Information” checkbox
- Select “Print” to print out the methods steps

The report can then be used to program the steps into the new firmware on the Profinia instrument through the touch screen interface.

Alternatively, or if you did not purchase Profinia software, you can record the method steps on Profinia planning worksheets found in Appendix K in the Profinia instruction manual.

For further assistance, contact Technical Support at 1-800-4BIORAD (1-800-424-6723) or lsg_techserv_us@bio-rad.com.

Windows is a trademark of Microsoft Corporation.

