

Autoimmune Panels

Ordering Information

Catalog No. Description

**BioPlex® 2200 ANA Screen with MDSS Panel**

665-1150	BioPlex® 2200 ANA Screen Reagent Pack .....	1 pack
663-1101	BioPlex® 2200 ANA Screen Calibrator Set .....	1 set
663-1131	BioPlex® 2200 ANA Screen Control Set .....	1 set

**BioPlex® 2200 Vasculitis Panel**

665-1850	BioPlex® 2200 Vasculitis Reagent Pack .....	1 pack
663-1800	BioPlex® 2200 Vasculitis Calibrator Set .....	1 set
663-1830	BioPlex® 2200 Vasculitis Control Set .....	1 set

**BioPlex® 2200 Anti-CCP**

665-3250	BioPlex® 2200 Anti-CCP Reagent Pack .....	1 pack
663-3200	BioPlex® 2200 Anti-CCP Calibrator Set .....	1 set
663-3230	BioPlex® 2200 Anti-CCP Control Set .....	1 set

**BioPlex® 2200 Antiphospholipid Syndrome Panels†**

665-1950	BioPlex® 2200 APLS IgG Reagent Pack .....	1 pack
663-1900	BioPlex® 2200 APLS IgG Calibrator Set .....	1 set
663-1930	BioPlex® 2200 APLS IgG Control Set .....	1 set
665-2050	BioPlex® 2200 APLS IgM Reagent Pack .....	1 pack
663-2000	BioPlex® 2200 APLS IgM Calibrator Set .....	1 set
663-2030	BioPlex® 2200 APLS IgM Control Set .....	1 set
665-2150	BioPlex® 2200 APLS IgA Reagent Pack .....	1 pack
663-2100	BioPlex® 2200 APLS IgA Calibrator Set .....	1 set
663-2130	BioPlex® 2200 APLS IgA Control Set .....	1 set

**BioPlex® 2200 Celiac Panels†**

665-2350	BioPlex® 2200 Celiac IgA Reagent Kit .....	1 pack
663-2300	BioPlex® 2200 Celiac IgA Calibrator Set .....	1 set
663-2330	BioPlex® 2200 Celiac IgA Control Set .....	1 set
665-2250	BioPlex® 2200 Celiac IgG Reagent Kit .....	1 pack
663-2200	BioPlex® 2200 Celiac IgG Calibrator Set .....	1 set
663-2230	BioPlex® 2200 Celiac IgG Control Set .....	1 set

† Available outside US only



# BioPlex® 2200 Autoimmune Panels

Bringing a new level of automation and productivity to autoimmune testing with a complete line of disease-specific multiplex panels

ANA Screen	APLS IgM†
Vasculitis	APLS IgA†
Anti-CCP	Celiac IgA†
APLS IgG†	Celiac IgG†

† Available outside US only




Bio-Rad Laboratories

For further information, please contact the Bio-Rad office nearest you or visit our website at [www.bio-rad.com/diagnostics](http://www.bio-rad.com/diagnostics)

Clinical Diagnostics Group

Website [www.bio-rad.com/diagnostics](http://www.bio-rad.com/diagnostics) U.S. 1-800-224-6723 Australia 61-2-9914-2800 Austria 43-1-877-8901 Belgium 32-9-385-5511 Brazil 5521-3237-9400 Canada 1-514-334-4372 China 86-21-61698500 Czech Republic 420-241-430-532 Denmark +45-4452-1000 Finland 358-9-804-22-00 France 33-1-47-95-60-00 Germany +49 (0)89-318-840 Greece 30-210-7774396 Hong Kong 852-2789-3300 Hungary +36-1-459-6100 India 1800-180-1224 Israel 972-3-9636050 Italy +39-02-216091 Japan 81-3-6361-7070 Korea 82-2-3473-4460 Mexico +52 (55)5488-7670 The Netherlands +31-318-540666 New Zealand 64-9-415-2280 Norway 47-23-38-41-30 Poland 48-22-3319999 Portugal 351-21-472-7700 Russia 7-495-721-14-04 Singapore 65-6415-3170 South Africa 27-11-442-85-08 Spain 34-91-590-5200 Sweden 46-8-555-127-00 Switzerland +41 (0)26-674-52-01/02/03 Taiwan 886-2-2578-7189 Thailand 662-651-8311 United Kingdom +44 (0)20-8328-2000





# BioPlex® 2200 Autoimmune Disease Panels

## Leading the way with an expanding multiplexed assay menu

Comprehensive assay menu offering a multiplexed solution for simultaneous detection and identification of the most clinically relevant autoantibodies.

### Offering the benefits of multiplex testing:

- High throughput multiplex testing improves workflow efficiencies and eliminates hands-on steps
- Each kit provides flexibility to test multiple autoantibodies using a small patient sample (5 µL)
- Consolidates traditional single-analyte tests into one panel
- Time-to-first result is approximately 45 minutes
- Ready-to-use multiplexed calibrators and controls, with curve stability up to 30 days

‡ Available outside US only

## Systemic Autoimmune Diseases

BioPlex® 2200 ANA Screen with MDSS Panel



The BioPlex® 2200 ANA Screen kit will help the clinician to better diagnose systemic autoimmune diseases by detecting only clinically significant autoantibodies.

### BioPlex® 2200 ANA Screen with MDSS

#### Analytes

- ANA Screen\*
- dsDNA (quantitative)
- Chromatin
- Ribosomal Protein
- SS-A (52 and 60)\*\*
- SS-B
- Sm
- Sm/RNP
- RNP (A and 68)\*\*
- Scl-70
- Jo-1
- Centromere B

\* ANA Screen is a composite roll up of all individual results  
\*\* Individual results available outside US

‡ Available outside US only



## Systemic Small Vessel Vasculitis

BioPlex® 2200 Vasculitis Panel



The BioPlex® 2200 Vasculitis kit will help the clinician diagnose some of the most severe autoimmune systemic vasculitic diseases affecting blood vessels in the lungs and kidneys.

### BioPlex® 2200 Vasculitis

#### Analytes

- Anti-MPO
- Anti-PR3
- Anti-GBM

## Antiphospholipid Syndrome

BioPlex® 2200 Antiphospholipid Syndrome (APLS) Panels†

### BioPlex® 2200 APLS IgG

#### Analytes

- Anti-Cardiolipin IgG
- Anti-Beta 2 Glycoprotein I (β2GPI) IgG

### BioPlex® 2200 APLS IgM

#### Analytes

- Anti-Cardiolipin IgM
- Anti-Beta 2 Glycoprotein I (β2GPI) IgM

### BioPlex® 2200 APLS IgA

#### Analytes

- Anti-Cardiolipin IgA
- Anti-Beta 2 Glycoprotein I (β2GPI) IgA

## Rheumatoid Arthritis

BioPlex® 2200 Anti-CCP



The BioPlex® 2200 Anti-CCP kit aids the clinician in the diagnosis of rheumatoid arthritis by providing important information regarding disease severity and prognosis.

### BioPlex® 2200 Anti-CCP

#### Analytes

- Anti-Cyclic Citrullinated Peptide (CCP) IgG

## Celiac Disease

BioPlex® 2200 Celiac Panels†

### BioPlex® 2200 Celiac IgA

#### Analytes

- Anti-tissue transglutaminase (tTG) IgA
- Anti-Deamidated Gliadin peptide IgA
- IgA Verification Bead, flags IgA-deficient samples

### BioPlex® 2200 Celiac IgG

#### Analytes

- Anti-tissue transglutaminase (tTG) IgG
- Anti-Deamidated Gliadin peptide IgG

