




KnowItAll U<sup>®</sup>

The World's Largest Collection of Spectra

**BIO-RAD**

# Campus-Wide Access to over 2 Million Spectra



"Can be utilized within the teaching curriculum, academic research in chemistry & biology, & in a library setting"

— Dr. Malcolm D' Souza,  
Professor of Chemistry,  
Wesley College (review in  
*Chemistry World*)

KnowItAll® U offers all faculty and students—at any computer—unlimited access to over 2 million spectra (with related properties and structures) including:

#### Bio-Rad Sadtler spectra

- 217,500 IR
- 68,900 <sup>1</sup>H NMR
- 52,900 <sup>13</sup>C NMR
- 2,700 MS
- 15,800 Raman
- 30,700 UV-Vis

#### John Wiley & Sons spectra

- 35,700 IR
- 3,800 NIR
- 162,000 <sup>1</sup>H NMR
- 203,000 <sup>13</sup>C NMR
- 42,000 XNMR
- 826,500 MS
- 6,400 Raman

#### Prof. Dr. Wolfgang Robien spectra (Universität Wien)


- 304,500 <sup>13</sup>CNMR
- 59,900 XNMR

#### Other

- 7,900 IR
- 2,400 Raman

**Unlimited simultaneous users** can search the KnowItAll U Spectral Library by full spectrum, peak, structure, substructure, property (name, molecular weight, CAS Registry Number, boiling point, etc.), or any combination.

### Search with KnowItAll AnyWare™ - A Simple Web Browser Interface



- Access the entire KnowItAll U data collection
- Simple, intuitive interface for users at any level
- Use with any operating system, including Windows, Macintosh, or Linux

## Search with KnowItAll Desktop Software For Advanced Analyses

### True Integration.


Instantly transfer data from one application to another.

### Wide Array of Chemistry Tools in One Interface.

Search, manage, and analyze spectroscopic, chromatographic, and chemical information; draw structures, create reports, and more!

### Intuitive Toolboxes.

Tools are logically grouped according to function/task.



### Licensing

KnowItAll U is provided as a site license on an annual subscription basis and is available to the faculty, staff, and students of degree-granting colleges and universities.

KnowItAll U may also include unlimited downloads of the KnowItAll desktop software throughout your institution, depending on the package chosen. In addition to spectral search tools, this software provides multiple tools for chemistry and related disciplines—all within a single user interface.

- Search KnowItAll U Data Collection
- Analyze & Process Spectra
- Spectral Identification
- Predict NMR Spectra
- Build & Manage User Databases
- Perform Chemometrics Analyses
- Draw Chemical Structures
- Create Reports



[www.knowitallu.com](http://www.knowitallu.com)

#### Supported Platforms / Operating Systems

KnowItAll AnyWare Web Browser

- Web browsers (including Google Chrome™, Internet Explorer®, Microsoft Edge®, Firefox®, and Safari™) on compatible operating systems (including Windows®, Macintosh®, and Linux®)

KnowItAll Informatics System Software

- Microsoft® Windows Vista, 7, 8, and 10
- Parallels Desktop for Macintosh



**Bio-Rad**  
Laboratories, Inc.