

Unity Interlaboratory Program

Optimize your laboratory performance as a member
of the world's largest community of QC users




Qualitative Reports

Blood Typing/Serum Indices Reports

Displays your laboratory data compared to your peer group data

- Supplies monthly and cumulative response distributions
- Compares your laboratory response distribution to the peer group distribution
- Displays visual histogram



Connectivity



Basic

WebConnect



Automated Uploads to Unity Software

A web-based connectivity solution that allows laboratories to easily upload QC data from LIS systems, middleware and/or instruments directly into Unity Real Time online or UnityWeb

- Eliminate manual keying of data
- Standard LIS QC reports can be used
- No software to install

Advanced

UnityConnect



A connectivity solution that allows QC data from LIS systems, middleware and/or instruments to be quickly and easily imported into Unity software and web services

- Eliminate manual keying of QC data
- Standard LIS QC reports can be used
- Automatic data import and real-time connection
- Transparent to the flow of data to the LIS
- Capture and analyze QC data economically from laboratory instruments not connected to the LIS
- Easy-to-install optional interface hardware
- Choose from soft, serial and ethernet connections

Ordering Information

Cat # Description

Unity Real Time LT

804	Installation Package	Single use
805-1	Unity Real Time LT	Annual subscription
825i	Unity Remote Installation & Training	Single use

UnityWeb

870-1	UnityWeb	Annual subscription
-------	--------------------	---------------------

Unity Real Time

804	Installation Package	Single use
804-1	Unity Real Time.	Annual subscription
825i	Unity Remote Installation & Training	Single use

Unity Real Time online

804-W1	Unity Real Time online	Annual subscription
--------	----------------------------------	---------------------

Bio-Rad Mission: Control

1100002	Mission: Control	Annual Subscription
1100005	Mission: Control	Subscription and Service Bundle

Cat # Description

Westgard Advisor

811-1	Westgard Advisor	Annual subscription
-------	----------------------------	---------------------

Unity Alert

806-1	Unity Alert.	Annual subscription
825i	Unity Remote Installation & Training	Single use

WebConnect

870-W1	WebConnect 2.0.	Annual subscription
825i	Unity Remote Installation & Training	Single use

UnityConnect

820-1	UnityConnect Software.	Annual subscription
825i	Unity Remote Installation & Training	Single use

Contact your Bio-Rad sales representative for region-specific information. Hardware solutions available for serial and network devices. Contact your Bio-Rad sales representative for recommendation.



Bio-Rad Laboratories, Inc.

For further information, please contact the Bio-Rad office nearest you or visit our website at www.bio-rad.com/qc-datamanagement

Clinical Diagnostics Group

Web site www.bio-rad.com/qualitycontrol **USA** 1 800 224 6723 **Australia** 61 2 9914 2800 **Austria** 43 1 877 8901 **Belgium** 32 03 710 53 00 **Brazil** 55 31 3689 6600 **Canada** 1 514 334 4372 **China** 86 21 61698500 **Czech Republic** 420 241 430 532 **Denmark** 45 4452 1000 **Finland** 358 9 804 22 00 **France** 33 1 47 95 60 00 **Germany** 49 0 89 318 840 **Greece** 30 210 7774396 **Hong Kong** 852 2789 3300 **Hungary** 36 1 459 6100 **India** 1800 180 1224 **Israel** 972 3 9636050 **Italy** 39 02 216091 **Japan** 81 3 6361 7070 **Korea** 82 2 3473 4460 **Mexico** 52 55 5488 7670 **The Netherlands** 31 318 540666 **New Zealand** 64 9 415 2280 **Norway** 47 23 38 41 30 **Poland** 48 22 3319999 **Portugal** 351 21 472-7700 **Russia** 7 495 721 1404 **Singapore** 65 6415 3170 **South Africa** 27 11 442 85 08 **Spain** 34 91 590 5200 **Sweden** 46 8 555 127 00 **Switzerland** 41 0 26 674 55 05 06 **Taiwan** 886 2 2578 7189 **Thailand** 662 651 8311 **United Kingdom** 44 0 20 8328 2000

Bio-Rad, Liquichek, Unity Real Time and UnityWeb are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. Acrobat is a trademark of Adobe Systems Incorporated.

© 2018 Bio-Rad Laboratories, Inc. Printed in USA QSD18-190 Q-1144 Rev. 07/2018

