



## TOX/SEE Rapid Urine Drug Screen Tests

When every minute counts, TOX/See is the answer.





# TOX/See Rapid Urine Drug Screen Test

Fast, flexible and convenient drug testing  
is more important than ever.

Results at  
5 minutes

TOX/See

vs.

Results in  
15 minutes

Competitor





## Walkaway Flexibility

Test results are stable and readable for one hour to accommodate busy laboratory schedules.

## Comprehensive Test Panels

Simultaneous detection of as many as 12 drugs per cassette.

## Rapid Results

Test results are available to read at 5 minutes.

## One-Step Test Protocol

One-step test protocol allows for easy training.

## Differentiates between Amphetamine, Methamphetamine and Ecstasy in One Step

Provides more specific information to the clinicians for better patient management and saves time and money for the laboratories.

## Specific Tests for Emerging Drugs of Abuse

**Buprenorphine** – Buprenorphine was recently approved as a substitution treatment for heroin addiction in the USA, following Europe, Australia and Singapore.<sup>1</sup> As the number of buprenorphine users increases, it will become more important to monitor treatment compliance and efficacy.

**Oxycodone** – From 2005 to 2009, oxycodone-related cases in the US emergency departments increased by 181%.<sup>2</sup> As this trend continues, the demand for prescription monitoring will increase.

1. Lee CE. (2006). Tackling Subutex abuse in Singapore. Singapore Med J 47(11), 919-21. <http://www.smj.org.sg/sites/default/files/4711/4711e1.pdf>. Accessed July 2, 2019.

2. Crane, E. H. (2015). The CBHSQ Report: Emergency Department Visits Involving Narcotic Pain Relievers. Rockville, MD: Substance Abuse and Mental Health Services Administration, Center for Behavioral Health Statistics and Quality. [https://www.samhsa.gov/data/sites/default/files/report\\_2083/ShortReport-2083.html](https://www.samhsa.gov/data/sites/default/files/report_2083/ShortReport-2083.html). Accessed July 2, 2019.

## Our Tests Feature a Combination of the Following Drugs at These Cutoffs

Amphetamine	AMP..... 1000 ng/mL	THC	THC ..... 50 ng/mL
Opiates	OPI ..... 300 ng/mL	Methamphetamine	MET ..... 1000 ng/mL
Opiates	OPI2000 ..... 2000 ng/mL	Barbiturates	BAR ..... 300 ng/mL
Benzodiazepines	BZO ..... 300 ng/mL	Tricyclic Antidepressants	TCA..... 1000 ng/mL
Methadone	MTD..... 300 ng/mL	Ecstasy	MDMA..... 500 ng/mL
Cocaine	COC ..... 300 ng/mL	Phencyclidine	PCP ..... 25 ng/mL
Oxycodone	OXY ..... 100 ng/mL	Propoxyphene	PPX..... 300 ng/mL
Buprenorphine	BUP ..... 10 ng/mL		

## Ordering Information

Catalog No.	Description	
194-5221	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MDMA, MET, MTD, OPI, OXY, PCP, TCA, THC.....	25 tests
194-5225	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MDMA, MET, MTD, OPI, PCP, THC.....	25 tests
194-5228	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MET, MTD, OPI, PCP, TCA, THC.....	25 tests
194-5230	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MET, OPI, PCP, TCA, THC.....	25 tests
194-5182	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MET, MTD, OPI, THC.....	25 tests
194-5183	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MTD, OPI 2000, PCP, THC.....	25 tests
194-5209	TOX/See Drug Screen Test AMP, BAR, BZO, COC, MTD, OPI, PCP, THC.....	25 tests
194-5216	TOX/See Drug Screen Test AMP, BZO, COC, OPI, THC.....	25 tests
194-5200	TOX/See Drug Screen Test AMP, COC, OPI, PCP, THC.....	25 tests
194-5184	TOX/See Drug Screen Test AMP.....	25 tests
194-5205	TOX/See Drug Screen Test BUP.....	25 tests
194-5188	TOX/See Drug Screen Test COC.....	25 tests
194-5224	TOX/See Drug Screen Test MDMA.....	25 tests
194-5190	TOX/See Drug Screen Test MET.....	25 tests
194-5191	TOX/See Drug Screen Test MTD.....	25 tests
194-5203	TOX/See Drug Screen Test OXY.....	25 tests
194-5194	TOX/See Drug Screen Test PCP.....	25 tests
194-5204	TOX/See Drug Screen Test PPX.....	25 tests
194-5234	TOX/See Drug Screen Test TCA.....	25 tests
194-5187	TOX/See Drug Screen Test THC.....	25 tests

BIO-RAD and TOX/SEE are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions.



**Bio-Rad  
Laboratories, Inc.**

For further information, please contact the Bio-Rad office nearest you  
or visit our website at [www.bio-rad.com/diagnostics](http://www.bio-rad.com/diagnostics)

Clinical  
Diagnostics Group

Website [www.bio-rad.com/diagnostics](http://www.bio-rad.com/diagnostics) U.S. 1-800-224-6723 Australia +61(2)9914-2800 Austria +43-1-877-8901 Belgium +32(3)710-53-00 Brazil +55(31)3689-6600  
Canada 1-514-334-4372 China +86 21 6169 8500 Czech Republic +420-241-430-532 Denmark +45-4452-1000 Finland +358-9-804-22-00 France +33-1-47-95-60-00  
Germany +49(0)89-318-840 Greece +30-210-7774396 Hong Kong +852-2789-3300 Hungary +36-1-459-6100 India +1800-180-1224 Israel +972-3-9636050 Italy +39-02-216091  
Japan +81-3-6361-7070 Korea +82-2-3473-4460 Mexico +52(55)5488-7670 The Netherlands +31-318-540666 New Zealand +64-9-415-2280 Norway +47-23-38-41-30  
Poland +48-22-3319999 Portugal +351-21-472-7700 Russia +7-495-721-1404 Singapore +65-6415-3170 South Africa +27-11-442-85-08 Spain +34-91-590-5200  
Sweden +46-8-555-127-00 Switzerland +41(0)26-674-55-05/06 Taiwan +886-2-2578-7189 Thailand +662-651-8311 United Kingdom +44(0)20-8328-2000