

Update to EN 61010-1 Amendment Statement of Operating Conditions

Catalog # Description

1652100 MicroPulser Electroporator

For research purposes only.

Specifications

See the following user manual updates according to the guidelines in IEC/EN/UL 61010-1/A1, section 5.4.2.

Environmental Conditions

Operating temperature 3.5–35°C, indoor use

Humidity 0–95%, in the absence of condensation

Altitude Up to 2,000 m

Pollution degree 2

Input Power

Supply voltage fluctuations ± 10% (unless otherwise specified),

Maximum current 2 A Maximum power 240 W

Overvoltage category II

BIO-RAD is a trademark of Bio-Rad Laboratories, Inc. All trademarks used herein are the property of their respective owner.